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---

 **SCHÜTZ DENTAL**  
Micerium Group

Schütz Dental GmbH · Dieselstr. 5-6 · 61191 Rosbach/Germany  
Telefon: +49 (0) 6003 814-0 · Telefax: +49 (0) 6003 814-906  
info@schuetz-dental.de · www.schuetz-dental.de

CE 0297

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If there is anything in this instruction for use that you do not understand, please contact our customer service department before using the product. As the manufacturer of this medical device, we inform our users and patients that all serious events occurring in connection with it must be reported to us (the manufacturers) as well as the relevant authorities in the Member State where the user and/or patient is resident.

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## Instructions for use

### **Purpose:**

Denture acrylics consist of a liquid component and a powder component. Dentures and splints are manufactured with these materials by the dental technician in the dental laboratory using various manufacturing processes. The denture acrylics are intended for the fabrication of splints, partial dentures, full dentures or model dentures.

### **General product description:**

Veins for denture acrylics on powder/liquid-basis for the individualization of dentures.

### **Composition:**

Cellulose fibres or polyamide, depending on the type

### **Indications:**

For aesthetic individualization of denture bases and for visualization of veins in the palate area.

### **Contra-indications / interactions:**

Do not use if there is a known allergy to any of the components. If a patient has known hyper-sensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. The dentist should consider known interactions and cross-reactions of the product with other materials already in the patient's mouth before using the product.

### **Application / processing:**

After weighing the liquid, add the desired amount of veins and mix homogeneously with a spatula. Then add the amount of acrylic powder weighed out according to the amount of liquid and proceed according to the processing instructions of the denture acrylic in question. The veins are suitable for both hot- and cold-curing denture base acrylics.

### **Storage:**

Storage temperature 10 - 25 °C / 50 - 77 °F.

### **Shelf-life:**

The maximum shelf-life is printed on the label of each packaging. Do not use after the expiration date.

### **Side-effects:**

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (e. g. allergies) or local discomfort, however, cannot be ruled out completely. All serious incidents occurring in connection with the use of this product must be reported to the manufacturer and the relevant competent authority.

### **Disposal:**

Veins can be disposed of as domestic waste. Residual quantities and packaging material must be disposed of in accordance with local and/or statutory regulations.