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# Alphalink Implant

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If there is anything in this instruction for use that you do not understand, please contact our customer service department before using the product.

As the manufacturer of this medical device, we inform our users and patients that all serious events occurring in connection with it must be reported to us (the manufacturers) as well as the relevant authorities in the Member State where the user and/or patient is resident.

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# Alphalink Implant

## Purpose

Fixing composites / root pin cements are self-curing or dual-curing two-component composites for dental use. They are used for the adhesive, permanent bonding of restorative materials to restorative materials, to the hard tooth structure or to the denture base, or to fix root pins inside the root canal.

## Product description

Alphalink Implant is a colorless purely chemically curing fixing composite for the adhesive connection of zirconium dioxide constructions on restorations made of zirconium dioxide, precious metals, non-precious metals and ceramics. The guidelines and specifications of DIN EN ISO 4049 apply.

## Designated users

For the application by dental technicians in a dental lab.

## Composition

- Monomer matrix: diurethane dimethacrylate, urethane acrylate
- Total filler content: 62% (by weight) anorganic fillers, glass filler (average particle size 5 µm), cristobalite powder (average particle size 3 µm), pyrogenic silica, dibenzoyl peroxide

## Indications

- permanently fixing indirect restorations made from metal (precious metal, non-precious metal, titanium), metal oxide ceramics (e. g. zirconium dioxide, aluminium oxide) or glass ceramics on implant abutments made from metal (e. g. titanium)
- bonding adhesive titanium bases to individual abutments made from zirconium dioxide

## Contraindications

- If a patient has known hyper-sensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision.
- Fixation of restorations on natural hard tooth substance.
- If not at least one of the components to be joined is made of zirconium dioxide.
- If a secure cleaning and drying process or the described application technique is impossible.
- Restorations which are to be fixed simultaneously on natural stumps and on implant posts

## Cross-reactions with other agents

- Phenolic substances (as e. g. eugenol) inhibit polymerization. For this reason, do not use liners containing such substances (e.g. zinc oxide eugenol cements).
- Oxidative disinfectants (e. g. hydrogen peroxide) may cross-react with the curing system and thus influence polymerization. Therefore, do not treat the automix cartridge with an oxidative disinfectant. Instead, disinfect by e. g. wiping with rubbing alcohol.
- When using this medical device, the attending dentist must consider commonly known cross-reactions or interactions of the medical device with other materials already present in the patient's mouth.

## A. Permanent bonding of indirect metal restorations (precious metal, non-precious metal, titanium), metal oxide ceramics (e. g. zirconium dioxide, aluminium oxide) or glass ceramics on implant abutments made from zirconium dioxide, zirconium dioxide constructions on implant abutments made from metal (e. g. titanium)

### A1. Removal of the temporary restoration

Remove the temporary restoration and, if necessary, remove the temporary abutment. Clean objects to be joined thoroughly with alcohol (e. g. isopropanol). Subsequently, dry with oil-free air. Avoid renewed moisture ingress and contamination.

### A2. Insertion and check of the permanent abutment

- Exclusively use an original center screw to fix the restoration permanently inside the patient's mouth. Adhere to the manufacturer's instructions.
- A retraction cord can be applied for better control of the restoration fit and subsequent excess removal of the fixing composite.

### A3. Fitting the restoration and drying the area of application

- Check the restoration's fit and occlusion.
- In case of brittle ceramic restorations, check the occlusion only very carefully as there is a danger of fracturing the unfixed object.
- If necessary, use a fine diamond tool at medium revolution and little pressure to carry out small corrections. Repolish any surfaces that have been ground.
- Relative drying of the working area – preferably with cotton rolls and parotid patches – is essential for adhesive cementation with composites. Absolute drying, e. g. by using a rubber dam, is recommended.

### A4. Pretreatment of the abutment surfaces

- Seal the screw canal of the abutment e.g. with temporary filling material to prevent the penetration of fixing material.
- Clean and dry the abutment (see A1.).
- Apply SEBOND Implant according to its instructions.

### A5. Pretreatment of the restoration surfaces

In general, follow the instructions of the framework material's manufacturer.

Other than that, the restorations are prepared as follows:

- Restorations made from metal, zirconium dioxide ceramics or aluminium oxide ceramics
  - Clean and dry the restoration (see A1.).
  - Sandblast the inside surfaces of the restoration with aluminium oxide (50 µm) at a maximum of 2 bar.
  - Dry the restoration with oil-free air.

**Please note: To obtain an ideal bond, do not clean the surfaces with phosphoric acid.**

Restorations made from lithium disilicate glass ceramics (e. g. IPS e.max Press, IPS e.max CAD)

- Etch with 5 % hydrofluoric acid for ca. 20 sec or in accordance with the material manufacturer's instructions.
- Rinse thoroughly with water and dry with an oil-free airflow.
- Apply SEBOND Implant in accordance with its instruction manual.

## A6. Application of Alphalink Implant

- Place a disposable mixing tip onto the dual-compartment cartridge. Discard the first bit of mixed material discharged from the tip.
- Press Alphalink Implant from the automix tip and apply the desired amount directly onto the restoration. Spread evenly.
- As the material will cure inside the used mixing tip, leave the tip on the cartridge as a closure until the next use.

## A7. Setting in the restoration and removal of any excess material

- Place the restoration in its destination with slight pressure. Fix it and hold it in place.
- Remove any excess fixing composite immediately with a disposable brush, cotton pellets, dental floss or an implant scaler. Take special care to remove excess material in difficult to reach areas (approximal, gingiva border) before the material cures.
- To avoid the formation of an oxygen inhibition layer, cover the margins of the restoration with glycerin gel/airblocker immediately after removing any excess material. This gel can be rinsed off with water after the fixing material has been cured completely.

## A8. Finishing the restoration

- If necessary, remove the retraction cord.
- Check the occlusion and the movements/function, correct if necessary.
- Polish all margins of the restoration with polishing tools or discs.

## B. Fixing individual zirconium dioxide abutments with adhesive titanium bases

### B1. Surface preparation of the adhesive titanium base

- Reduce the flare of the titanium base to the implant diameter in the area of the bonding (transition zirconium dioxide / titanium) according to subgingival anatomical aspects. Take into account the minimum wall thickness of the zirconium dioxide portion.
- Protect the hexagonal connection of the titanium base by screwing it onto a lab implant for the time of the bonding procedure.

### Important information: Only use a lab screw for fixation!

The scan screw is not suited to be used in the bonding procedure. If the bonding material penetrates the microgap, the screw cannot be removed anylonger.

- Block out the screw thread in the titanium base with wax.
- Clean the titanium bases thoroughly with alcohol (e. g. isopropanol) and dry with an oil-free airflow. Avoid renewed moisture ingress and contamination.
- Carefully sandblast the adhesive titanium base with 50 µm aluminium oxide at a maximum of 2.5 bar.
- Apply SEBOND Implant according to the manufacturer's instructions.

### B2. Surface treatment of the zirconium dioxide abutment

- Clean the abutments thoroughly with alcohol (e. g. isopropanol) and dry with an oil-free airflow. Avoid renewed moisture ingress and contamination.
- Carefully sandblast the inside of the abutment with 50 µm aluminium oxide at a maximum of 1.0 bar.
- Apply SEBOND Implant according to the manufacturer's instructions.

## B3. Bonding the abutment and removal of any excess material

- Place a disposable automix tip onto the dual-compartment cartridge. Discard the first bit of mixed material discharged from the tip.
- Press Alphalink Implant from the automix tip and apply the desired amount directly onto the abutment. Spread evenly.
- As the material will cure inside the used mixing tip, leave the tip on the cartridge as a closure until the next use.
- With help of a microscope, carefully and thoroughly remove any excess adhesive from the screw channel.
- If necessary, use a suitable reamer to smoothen the inside of the screw channel in the titanium insert.

## B4. Finishing the fixed abutment

- After bonding, the tapered zirconium dioxide abutments are reworked in the milling unit under water cooling until a silky-matt gloss is achieved. Finally, the abutments are polished to a high gloss by hand.

## Troubleshooting / FAQ

Error	Cause	Remedy
Bond does not last after use of cement	Bonding surface contaminated	Ensure contamination-free work after sandblasting and primer application
	Components of the cement unevenly mixed (air bubbles)	Use mixing tip and discard the first ca. 0.25 g discharged from the tip
	Processing time exceeded during application	Observe processing time and mix only small quantities
Cement does not cure	Use of self-etch bonding systems	DDo not use any surface primers foreign to the system. Only use original primers.
	Contamination of the bonding surfaces with materials containing eugenol	Use only eugenol-free products.
	Incorrect storage of the cement (temperature too high)	New purchase of cement and storage in a refridgerator.
	Mixing of uneven quantities of the two cement components	Use original mixing tip and discard the first bit of material discharged.

## Processing time

The processing and curing times depend on the ambient temperature. As soon as Alphalink Implant is discharged from the automix tip, the following times apply:

	Room temperature ca. 21 °C (ca. 69 °F)	Intraoral temperature 37 °C (98.6 °F)
Processing time	2–3 min	ca. 1 min
Curing time, incl. processing time	ca. 5 min	2–3 min

## Contraindications / interactions

If a patient has known hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. Before using the product, the dentist should consider known interactions and cross-reactions of the product with other materials already in the patient's mouth.

Unpolymerized acrylics may lead to allergic reactions of the skin. Therefore, the user should take suitable precautions. In case of irritation or known allergy to any of the substances listed in the composition, do not use.

## Mixing ratio

By using the automix cartridge with the mixing tip, Alphalink Implant is always mixed in the ideal ratio of 1 : 1.

## Storage and handling

- Alphalink Implant must be refridgerated. We recommend a storage temperature of 3 °C – 9 °C (ca. 37 °F – 48 °F). Keep the used automix tip on the cartridge for storage.
- The paste should reach room temperature before processing, as refrigerator temperature may make pressing out and mixing more difficult. Alphalink Implant should be processed quickly after discharging from the automix cartridge and the restoration should be placed quickly.

## Shelf-life

The maximum shelf-life is printed on the label of each primary packaging. Do not use after the expiration date.

## Note on side effects

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (e. g. allergy) or local discomfort, however, cannot be ruled out completely. All serious incidents arising in connection with the use of this product must be reported to the manufacturer indicated below and to the relevant competent authority.

## Disposal

Residual quantities and packaging materials are to be disposed of in accordance to local and/or legal requirements.

CE 0297

**SCHÜTZ DENTAL**  
Micerium Group

Schütz Dental GmbH  
Dieselstr. 5-6 · 61191 Rosbach/Germany  
Tel.: +49 (0) 6003 814-0 · Fax: +49 (0) 6003 814-906  
info@schuetz-dental.de · www.schuetz-dental.de