# Capo Bulk Fill Composite

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### Instructions for use

Capo Bulk Fill Composite is a light curing, flowable, highly radiopaque (210% A) composite for posterior restorations, for direct filling and the „bulk fill“ technique. Capo Bulk Fill Composite is available in syringes and compules. The compules are for single use. Please do not reuse them, as this makes it impossible to rule out contamination and germ formation. The product complies with the guideline DIN EN ISO 4049.

**Composition:**  
**Monomer matrix:**  
Diurethane dimethacrylate, 1,4-butanediol dimethacrylate  
**Total filler:**  
77 % by weight (57 % by volume) inorganic filler (0.04-15 µm)

#### Indications:

- Fillings with layer thicknesses to 4 mm in Black’s class I cavities.
- Base fillings/liners with layer thicknesses up to 4 mm in Black’s class II cavities. For the following layers in a class II cavity, only use composites which are suited to filling Black’s class II cavities.
- Cavity lining – as a first (bottom) layer in Black’s class I and II cavities.

#### Application:

##### Preparatory measures

Before beginning the treatment, clean the hard tooth substance with non-fluoride polishing paste.

##### 1. Cavity preparation

Prepare the cavity with a minimally invasive method as generally required for adhesive techniques. All enamel margins in the anterior region must be bevelled. Do not bevel the margins in the posterior region and avoid slice preparations. Spray the cavity with water to clean it, remove all debris and dry it. Isolate the cavity to keep it dry. It is advisable to place a rubber dam.

##### 2. Approximal contact areas

For cavities with approximal sections, place a transparent matrix and fix it in place.

##### 3. Adhesive system

Etch (e.g. Capo Etch) and bond (e.g. Capo Bond) according to manufacturer’s instructions. Avoid drying out the dentin. Dried, etched surfaces appear chalky-white and must not be contaminated before the application of a bonding agent. In case of contamination with saliva, rinse and dry again. If necessary, etch once more.

##### 4a. Application of Capo Bulk Fill Composite (syringes)

Lock the application tip on the syringe by turning it clock-wise. Place layers of max. 4 mm of Capo Bulk Fill Composite directly into the cavity using the curved application tips supplied with the material. For hygiene reasons, these application tips must only be used once and then discarded. Avoid air bubbles when applying the material. Wet the prepared tooth surfaces thoroughly. Polymerize each layer with a conventional polymerization unit for 40 seconds, holding the light guide as close as possible to the surface. Check the compatibility of the polymerization unit. Ensure that the polymerization time is sufficient. Test the light output before each use. During polymerization, a dispersion layer forms on the material surface. This layer must not be touched or removed. It forms the bond with the layer of composite applied subsequently (e.g. with a universal or posterior composite) or for creating aesthetic cusps. Alternatively, Capo Bulk Fill Composite can be applied in a thin layer to form a liner on exposed dentin. The material can be polymerized in layers of up to 4 mm thickness.

##### 4b. Application of Capo Bulk Fill Composite (compules)

Insert the compule into the dispenser. Remove the sealing cap. Place the compule in the correct angle towards the cavity. Inject the material into the cavity. Apply slow and steady pressure to the compule. Do not use excessive force! To remove the compule from the dispenser after use, retract the plunger. Next, remove the compule.  
**Please note:** For hygiene reasons, compules are intended for single use only.

##### 5. Finishing

Capo Bulk Fill Composite can be finished and polished directly after polymerization. To finish, use finishing diamond tools, flexible discs, silicone polishers and polishing brushes. Check occlusion and articulation. Grind down if necessary, thus avoiding premature or defective contacts on the surface of the filling.

#### Please note

- When creating large and time-consuming restorations either cover the material with a lightproof foil or temporarily move the operating light away from the site in order to avoid premature polymerization of the material.
- For hygienic reasons, the curved application tips supplied with the material must only be used once, as a contamination of the material and a transfer of germs cannot be ruled out.
- Use a polymerization unit with an emission spectrum of 350–500 nm to cure this material. As the required physical properties can only be achieved if the lamp is functioning correctly, its luminous intensity must be checked regularly as described by the lamp’s manufacturer.

#### Hazard statement

May cause an allergic skin reaction.

#### Precautionary statements

Wear protective gloves/protective clothing/eye protection/face protection. If skin irritation or rash occurs: Get medical advice/ attention.

#### Contraindications:

Unpolymerized composite may cause skin allergies. The user must take adequate precautions (e.g. protective gloves). Do not use this material, if the user or patient should have known allergies against one of the components listed under „composition“.

#### Side-effects:

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us. To prevent possible reactions of the pulp in cavities where the dentin is exposed, the pulp must be protected adequately (apply e.g. a calcium hydroxide preparation).

#### Interaction with other substances:

As phenolic substances (such as eugenol) inhibit polymerization, do not use cavity liners containing such substances (e.g. zinc-eugenol-cements).

#### Desinfection/Protection from cross-contamination

Place the functioning syringe with attached application tip into a suitably shaped barrier sheath; pierce end of sheath with cannula, exposing the cannula for use. Using a barrier sheath facilitates cleaning and disinfection of the syringe between patients. After use of sheathed syringe, remove delivery tip and sheath by grasping on the hub of the delivery tip through the sheath, twist and remove tip along with sheath. Discard used tip and sheath in appropriate waste stream. Replace syringe storage cap.

**Disinfect:** After removal of the delivery tip and sheath, disinfect this product using an intermediate level disinfection process (liquid contact) as recommended by the Centers for Disease Control and endorsed by the American Dental Association, Guidelines for Infection Control in Dental Healtho-Care Settings – 2003 (Vol. 52; No. RR-17), Centers for Disease Control and Prevention (USA).

#### Storage:

Store at 10–25°C (50–77°F). Avoid direct sunlight. Close syringes tightly immediately after use. Let the material reach room temperature before use. After use, retract the plunger of the syringe slightly to prevent the apertures from becoming blocked. Do not use after the expiry date (refer to label on syringe). For use by dentists only. Keep

out of children’s reach. This product was developed specifically for the described range of application.

It must only be used according to the above instructions. The manufacturer will not be held liable for damages caused by in-correct handling or incorrect processing of the material.

#### Contains 1,4-butanediol dimethacrylate

**Warning:** May cause an allergic skin reaction. Avoid breathing vapours/spray. Wear protective gloves/protective clothing/eye protection/face protection. If skin irritation or rash occurs: Get medical advice/attention.

#### Troubleshooting

Trouble	Cause	Remedy
composite does not cure	luminous intensity of the polymerization unit insufficient	check luminous intensity; replace light source, if necessary
	emitted spectral range of the polymerization unit insufficient	consult manufacturer of polymerization unit; recommended spectral range: 350 – 500 nm
composite seems to be too hard/firm inside the syringe	material was stored at temperatures below 10°C for a longer period of time	let composite reach room temperature before use
	syringe was not closed tightly which caused part of the material to cure	close syringe tightly after each use
composite does not cure sufficiently	layer thickness per polymerization cycle too high	keep to max. layer thickness of 4 mm
restoration seems too yellow when compared to color reference	insufficient polymerization of the composite layers	repeat polymerization cycle several times, for a minimum of 40 sec.

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