CapoFlow

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If there is anything in this instruction for use that you do not understand, please contact our customer service department before using the product.
As the manufacturer of this medical device, we inform our users and patients that all serious events occurring in connection with it must be reported to us (the manufacturers) as well as the relevant authorities in the Member State where the user and/or patient is resident.

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Instructions CapoFlow

CapoFlow is a light curing, flowable, radioopaque, low viscosity composite. Complies with the guidelines described in DIN EN ISO 4049.

Composition

Glass powder, di-urethane dimethacrylate, Bis-GMA, tetramethylene dimethacrylate, splitter polymer

Total filler:

55% by weight (36% by volume) inorganic fillers (0.005 – 3 µm)

Indications

- Fissure sealing
- Extended fissure sealing on molars and premolars
- Fillings in Black's class V cavities (cervical caries, eroded areas in roots, wedge-shaped defects)
- Minimally invasive fillings in Black's class I and II cavities in areas not exposed to severe occlusal loads
- Minimally invasive fillings in Black's class III cavities
- Restoring defects in enamel
- Blocking out undercut
- Minimal adjustments to the contours and shade of the enamel

Application

Preparatory measures

Before commencing the treatment, clean the tooth with non-fluoride polishing paste. Use a Vita® shade guide to select the shade while the tooth is still moist.

1. Cavity preparation

Prepare the cavity minimally invasively as generally required for adhesive techniques. All enamel margins in the anterior region must be bevelled. Do not bevel the margins in the posterior region and avoid slice preparations. Spray the cavity with water to clean ist, remove all debris and dry it. The cavity must be isolated. It is advisable to place a rubber dam.

2. Pulp protection /Cavity liner

If an enamel-dentin adhesive is used, no cavity liner is required. In very deep cavities those areas in close proximity to the pulp must be coated with a calcium hydroxide material.

3. Approximal contact areas

For cavities with approximal sections, place a transparent matrix and fix it in place.

4. Adhesive system

Etch (e. g. with Capo Etch) and bond (e. g. with Capo Bond) according to the manufacturer's instructions.

5. Applying CapoFlow

Place thin layers (max. 2 mm) of CapoFlow directly in the cavity using the curved application tips supplied with the material – for hygienic reasons, they are for single-use only. When applying the material with the tip, ensure that no air bubbles become entrapped. Ensure that the prepared tooth surfaces are wetted thoroughly. Light cure each layer with a commercially available polymerization unit for 40 seconds, holding the light guide as close as possible to the surface. During polymerization, a dispersion layer forms on the material surface. This layer must not be touched or removed. If forms the bond with the layer of composite applied later.

6. Trimming

CapoFlow can be trimmed and polished immediately after curing using finishing diamonds, flexible disks, silicone polishers and polishing brushes.

Please note

- When placing time-consuming restorations, move away the dental light temporarily or cover the material with an opaque foil to prevent premature polymerization of the composite.
- For hygienic reasons, the curved application tips supplied with the material must only be used once!
- Use a light polymerization system with an emission range of 350 – 500 nm to polymerize the material. The required physical properties are only reached if the polymerization light functions properly. Therefore, it is necessary to check the light intensity regularly according to the manufacturer's instructions.

Light intensity for polymerization	≥ 650 mW/cm²
Wave length for polymerization	350–500 nm
Polymerization time	40 sec.

Hazard statement / precautionary statement
Contains tetramethylene dimethacrylate, di-urethane dimethacrylate
Warning: May cause an allergic skin reaction. Harmful to aquatic life with long lasting effects.
Avoid breathing vapours/spray. Avoid release to the environment. Wear protective gloves. If skin irritation or rash occurs: Get medical advice/attention.

Contra-indications

If a patient has known hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. In such cases, we will supply the composition of our medical device upon request. The dentist should consider known interactions and cross reactions of the product with other materials already in the patient's mouth before using the product.

Side-effects

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us.

In case of exposed dentin in a cavity, ensure a proper protection of the pulp to avoid possible pulp reactions (e. g. line the cavity with a calcium hydroxide compound).

Interaction with other substances

As phenolic substances (such as eugenol) inhibit polymerization, do not use cavity liners containing such substances (e. g. zinc-eugenol-cements).

Disinfection / Protection from cross-contamination:

Place the functioning syringe with attached application tip into a suitably shaped barrier sheath; pierce end of sheath with cannula, exposing the annula for use. Using a barrier sheath facilitates cleaning and disinfection of the syringe between patients.

After use of sheathed syringe, remove delivery tip and sheath by grasping on the hub of the delivery tip through the sheath, twist and remove tip along with sheath. Discard used tip and sheath in appropriate waste stream. Replace syringe storage cap.

Disinfect – After removal of the delivery tip and sheath, disinfect this product using an intermediate level disinfection process (liquid contact) as recommended by the Centers for Disease Control and endorsed by the American Dental Association, Guidelines for Infection Control in Dental Health-Care Settings – 2003 (Vol. 52; No. RR-17), Centers for Disease Control and Prevention (USA).

Storage

Store at temperatures between 10 °C and 25 °C (50 °F to 77 °F). Screw the cap back onto the syringe tightly after each use. Let the material reach room temperature before use. Withdraw the plunger slightly after use to keep the outlet from becoming plugged. Do not use after the expiration date (see label on the syringe). Only for use in dentistry. Keep out of children's reach. This product has been developed for the specific use illustrated above. Only process as described in these instructions. The manufacturer will not be held liable for any damages that result from improper use or improper processing.

Troubleshooting

Trouble	Cause	Remedy
Composite does not cure	Luminous intensity of the polymerization unit insufficient	Check luminous intensity; replace light source, if necessary
	Emitted spectral range of the polymerization unit insufficient	Consult manufacturer of polymerization unit; recommended spectral range: 350 – 500 nm
Composite seems to be too hard/ firm inside the syringe	Material was stored at temperatures below 10°C for a longer period of time	Let composite reach room temperature before use
	Syringe was not closed tightly which caused part of the material to cure	Close syringe correctly with the cap after each use
Composite does not cure sufficiently	Layer thickness per polymerization cycle too high	Keep to max. layer thickness of 2 mm
Restoration seems too yellow when compared to color reference	Insufficient polymerization of the composite layers	Repeat polymerization cycle several times, for a minimum of 40 sec.

*Vita is a registered trademark of Vita® Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen, Germany.

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*Vita è un marchio registrato della Vita® Zahnfabrik H. Rauter GmbH & Co. KG., Bad Säckingen