Capo Hybrid

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Capo Hybrid Instructions for use

Capo Hybrid is a light curing, hybrid composite containing an ultrafine, radiopaque glass filler and is indicated for placing fillings using adhesive techniques. It can be polished to a high lustre. Due to the ultrafine particle filler, extremely homogeneous restorations can be placed which are easily polished to a high lustre. The chameleon effect matches the shade of the filling perfectly to the tooth structure. The guidelines of EN ISO 4049 have been complied with. Capo Hybrid is supplied in practical 5 g screw syringes or **single-use**, 0.3 g “Fillis”. Do not use fillis more than once. Otherwise, a contamination of the material or a possible transfer of germs cannot be excluded. Capo Hybrid should be used with Capo Etch gel and Capo Bond which belong to the same system.

Intended use

Light-curing material for filling dental cavities.

Composition
Monomer matrix
Glass powder, aliphatic urethane dimethacrylate, silicon dioxide, Bis-GMA, tetramethylene dimethacrylate
Total filler
75% by weight (58% by volume) anorganic filler (0.005 – 3.0 µm)

Indications

- Direct anterior and posterior restorations in Black’s class I, II, III, IV, and V cavities.
- Indirect restorations such as inlays, onlays and laminate veneers
- Extended fissure sealing in molars and premolars
- Endodontic posts
- Splinting mobile teeth
- Adjusting the contours and shades to improve aesthetics

Application

Pre-treatment

Before commencing the treatment, clean the tooth with non-fluoride polishing paste. Use a Vita® shade guide to select the shade while the tooth is still moist.

- Cavity preparation**
Minimally-invasive preparation of the cavity as generally required for adhesive techniques. All enamel margins in the anterior region must be bevelled. Do not bevel the margins in the posterior region and avoid slice preparations. Spray the cavity with water to clean it, remove all debris and dry it. The cavity must be isolated. It is advisable to place a rubber dam.
- Pulp protection / Cavity liner**
If an enamel-dentine adhesive is used, no cavity liner is required. In very deep cavities those areas in close proximity to the pulp must be coated with a calcium hydroxide material.
- Approximal contact areas**
When filling cavities with approximal sections, place a transparent matrix and fix it in place.
- Adhesive system**
Ätch (e.g. with Capo Etch) and bond (e.g. with Capo Bond) in accordance to the manufacturer’s instructions.
- Placing the composite in the cavity**
 - Compulas**
Insert the compules in the dispenser and remove the sealing cap. Exert pressure slowly and uniformly to dispense the required amount of filling material from the compules directly into the cavity, layer-by-layer. Contour the composite with standard metal instruments. Do not exceed a layer thickness of 2 mm.
Please note: For hygienic reasons compules are for single-use only.

b) Screw syringes

Extrude the required amount of composite from the screw syringe, place it in the cavity with a standard metal instrument and contour it. The layers must not be thicker than 2 mm.

6. Polymerization

Each layer must be polymerized for 40 sec. in a conventional halogen polymerization system or in an LED polymerization system. As an alternative, polymerize the material 2 x 3 sec. with a plasma polymerization system. Hold the light guide as close to the surface of the filling as possible. In case of fillings with several surfaces, polymerize each surface separately. Due to the effect of the oxygen in the air, a thin smear layer of unpolymerized material remains on the surface of each layer. This bonds the layers chemically and must not be touched or contaminated with moisture.

7. Trimming

Capo Hybrid can be trimmed and polished immediately after polymerization using finishing diamonds, flexible discs, silicone polishers and polishing brushes. Check the occlusion and articulation and spot grind to eliminate high spots or undesirable paths of articulation from the surface of the filling.

Indirect method

Cavity preparation

The cavity should be prepared as minimally invasively as possible with only slightly diverging sides. To prevent the material fracturing, the layer must have a minimum thickness of 1.5 mm in the lateral and vertical aspects. All internal edges and angles must be rounded. Avoid slice preparations. Prepare a flat cervical shoulder – do not bevel it. Any un-avoidable undercuts must be blocked out with glass ionomer cement. Use slightly tapered diamonds with rounded tips for the preparation. Coat those areas of dentine in close proximity to the pulp with a thin layer of calcium hydro-xide material. Cavity liners containing eugenol are contraindicated.

Impression and temporary restoration

Once the impression has been taken, a composite temporary restoration is fabricated (e.g. with Temdent, Schütz Dental). This may only be cemented with a non-eugenol cement.

Fabrication of inlays, onlays and veneers

Cast the impression with hard stone plaster. Allow the model to set and pull off the impression. Block out the undercuts and apply an oil-free separating agent to the model. Build up the inlay on the model layer-by-layer. Build up the approximal and deep occlusal sections first. The layer thickness of each individual layer must not exceed 2 mm. Each layer is polymerized separately with a conventional polymerization lamp (e.g. HiLite Power, Heraeus Kulzer, intermediate polymerization 90 sec., final polymerization 180 sec.). The restoration is lifted from the model. It is finished and polished to a high shine. Clean the restoration thoroughly with soap and water. Rinse with air/water spray and dry.

Insertion of the inlay, onlay or veneer

Remove the temporary restoration and clean the cavity. Place a rubber dam before cleaning and drying the prepared surfaces of the tooth. Exert only gentle pressure on the restoration to check its fit. Do not use excessive force. If necessary, trim those surfaces facing the cavity to improve the fit. To avoid fractures, do not check the occlusion during the fitting.

Etch (e.g. Capo Etch) and bond (e.g. Capo Bond) according to the manufacturer’s instructions.

Fixing the restoration

The restoration is fixed with a conventional, dual-curing fixing composite (e.g. Alphalink Cem, Schütz Dental). Please adhere to the manufacturer’s instructions.

Special notes

- The working time under a surgical lamp is approximately 2 minutes.
- In case of time-consuming restorations, the surgical lamp should be either temporarily moved away from the working area or the material should be covered by an opaque foil in order to prevent the composite from curing too early.
- Use a light-curing unit with an emission spectrum of 350 – 500 nm for the polymerization of this material. As the required physical properties can only be achieved when the lamp works correctly, its luminous intensity must be checked regularly as described by the manufacturer.

Light intensity for polymerization	650 mW/cm²
Wavelength for polymerization	350–500 nm
Polymerization time	40 sec.

Warnings and safety instructions

Contains tetramethylen dimethacrylate, aliphatic urethane dimethacrylate

Caution: May cause allergic skin reactions. Harmful to aquatic organisms, with long lasting effects. Avoid inhalation of vapor /aerosol. Wear protective gloves. In case of skin irritation or rash: Get medical advice/ attention.

Contraindications / interactions

If a patient has know allergies against or hypersensivities to a component of this product, it may not be used or only under strict medical supervision by the doctor/dentist. The dentist should consider known interactions and cross-reactions of the product with other materials already in the patient’s mouth before using the product.

Side-effects

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us.

To prevent possible reactions of the pulp in cavities where the dentine is exposed, the pulp must be protected adequately (e.g. calcium hydroxide preparation).

Interactions with other substances

As phenolic substances (such as eugenol) inhibit polymerization, do not use cavity liners (such as zinc-oxide eugenol cements) containing such substances.

Storage

Store at 10–25 °C (50–77 °F).

Close the screw syringes tightly immediately after use. The material should be at room temperature before use. Retract the plunger of the syringe slightly to prevent the apertures becoming blocked.

Do not use after expiry date (refer to label on syringe). For use by dentists only. Keep out of reach of children. This product was developed specifically for the described range of applications. It must be used as described in the instructions. The manufacturer is not liable for damage caused by handling or processing the material incorrectly.

Package contents	
Capo Hybrid Set (Fillis):	Art. No.: 232110
6 x 10 Fillis, 0.3 g each	

Capo Hybrid Refills (Fillis):	
20 Fillis x 0.3 g each	

Capo Hybrid Set (screw syringes):	Art. No.: 232011
6 screw syringes	

Capo Hybrid Refills (screw syringes):	
5 g screw syringe	

* Vita is a registered trademark of Vita® Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen, Germany.

Problem	Cause	Remedy
Composite does not polymerize properly	Light output of the polymerization lamp is inadequate <p>Emitted wavelength range of the polymerization lamp is inadequate</p>	Check the light output and change the light source if required <p>Consult the manufacturer of the polymerization lamp. Recommended wavelength range: 350 – 500 nm</p>
Composite in the syringe is sticky and soft, colorless liquid separates in the syringe	Material has been stored for a longer period ≥ 25 °C <p>Material has been kept in a syringe warmer for too long</p>	Observe the storage temperature. Store at 10 – 25 °C. <p>Do not keep syringes in a syringe warmer for longer than one hour per session.</p>
Composite appears too hard and firm in the syringe	Material was stored at < 10 °C for a longer period of time. <p>Syringe not properly sealed, composite partially cured</p>	Allow the composite to reach room temperature before use; use a syringe warmer if necessary <p>Always seal the syringe properly with the cap after taking out composite</p>
Inlay /onlay is not properly retained when fitted	Restoration is too opaque to be cemented using only light-curing composite	Use dual-curing luting composite
Composite does not cure completely (dark or opaque shades)	Composite layers applied too thickly for each polymerization cycle	Adhere to a max. thickness of 2.0 mm per layer
Restoration appears too yellow compared with the shade guide	Inadequate polymerization of the composite layer	Repeat the exposure cycle several times; min. 40 seconds.

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* Vita è un marchio registrato della Vita® Zahnfabrik H. Rauter GmbH & Co. KG., Bad Säckingen