Capo Bond Activator

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If there is anything in this instruction for use that you do not understand, please contact our customer service department before using the product. As the manufacturer of this medical device, we inform our users and patients that all serious events occurring in connection with it must be reported to us (the manufacturers) as well as the relevant authorities in the Member State where the user and/or patient is resident.

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Instructions for use

General product description:
Capo Bond Activator is an additive which is simply mixed with Capo Bond to convert it to a dual-curing or self-curing primer and bonding system.
Capo Bond Activator is highly recommended when using a dual-cure or self-cure primer and bonding system.

Composition:
Ethanol, methacrylate carboxylic acid ester, initiators

Indications:

- As an additive to Capo Bond for the strong bond of restorative composites on natural hard tooth substances, such as enamel and dentin.
- As an additive to Capo Bond for the fixation of indirect restorations, e.g. ceramic inlays and composite inlays, onlays, veneers, crowns and bridges, on natural tooth structure with a chemically curing or dual-curing luting composite.

Contraindications:

- Do not use in case of a known allergy against any of the components.
- Unpolymerized bonding may cause allergic reactions of the skin. Users should wear protective gloves.

Application / processing of Capo Bond Activator as an additive to Capo Bond:
Processing of Capo Bond Activator as an additive to Capo Bond for the firm bonding of restorative composites to natural tooth structure and for the fixation of indirect restorations to natural tooth structure.

Preparatory tasks:

Insulation and preparation of the cavity:
The use of a rubber dam is recommended. Before starting the preparation, clean the tooth to be treated with pumice powder and water. Prepare the cavity with a minimal tooth reduction to avoid losing too much healthy tooth substance. The cavity margin must be beveled slightly (0.5 - 1.0 mm) at the enamel to increase the bonding surface of the enamel and thus the bond strength.

Protection of the pulp:
In case of deep cavities that come close to the pulp, cover the according surfaces with a thin layer of a curing calcium hydroxide compound.

Etching the hard tooth substance:
The total-etch technique is recommended: Apply a phosphoric acid based etching gel, e.g. Capo Etch, to the enamel and dentin surfaces according to the manufacturer's instructions. Do not let the dentin dry out, because a slightly moist dentin surface is important for the effectiveness of the mixture of Capo Bond and Capo Bond Activator. After the treatment, the etched enamel surface appears distinctly chalky-white. The etched hard tooth substance must neither be touched nor contaminated with saliva before applying the mixture. In case of any contamination, the enamel surface must be etched again, flushed with water and dried.

Preparation and application of the mixture of Capo Bond and Capo Bond Activator:

Before use, shake Capo Bond well. Drip one drop of Capo Bond onto a mixing plate. Next, add one drop of Capo Bond Activator. Mix for 15 seconds under dimmed light.

Application on dentin and enamel:
Use a brush to apply a sufficient quantity of the mixture of Capo Bond and Capo Bond Activator onto the dentin surfaces and enamel surfaces. Brush intensively for 30 seconds to create a homogenous layer. Ensure that a uniform layer has been applied and that the treated dentin and enamel surfaces are kept moist with the mixture for the prescribed time. Remove any surplus material. Remove all volatile components by gently blowing with compressed air (approx. 15 seconds), spreading the adhesive to form an even layer. Do not dry out the dentin in the process.

Optionally, the wetted surfaces can be exposed to a suitable polymerization lamp for 20 seconds before a second layer is applied. If the layer is not exposed to a polymerization lamp, it must be allowed to cure for 3 minutes.

Application of a second layer:
Once again, apply a sufficient quantity of the dual-curing mixture of Capo Bond and Capo Bond Activator to the surfaces to be bonded. Brush intensively for 30 seconds. Next, dry with oil-free compressed air for ca. 15 seconds. The dual-curing system will polymerize on its own within 3 minutes afters its application and and the intermediate insertion of the restoration.

Placing the filling, polymerization and finishing of filling composites, application of the luting composite:

After the application of the dual-curing mixture of Capo Bond and Capo Bond Activator, the light-curing composite filling can be placed or the filling can be fixed with a light-curing or dual-curing luting composite. For the placing of the filling (e. g. with Capo Hybrid Composite), for the polymerization and for finishing the composite material, please adhere to the manufacturer's instructions. Reference is also made to the manufacturer's

instructions for use for the application and curing of luting composites.

Error	Cause	Remedy
mixture of Capo Bond and Capo Bond Activator does not cure	Light output of the polymerization lamp not sufficient <p>mixing ratio is not 1:1</p>	check the light output; replace the light source, if necessary <p>observe the mixing ratio</p>
restoration does not stay in place	working times and especially etching times were not observed	observe the working and etching times quoted
	individual layers were not polymerized separately	separate polymerization of each layer
	mixture was not brushed intensively	carefully observe the application instructions
material dried out or cured	unsuitable storage conditions, bottle not closed correctly	screw the lid back onto the bottle after each time material is removed
no homo-genous layer	mixture applied irregularly and with a non-sufficient layer thickness	apply fresh material

Guarantee:
Our processing instructions, regardless of whether they are issued verbally, in writing or in the form of a practical demonstration, are based upon our own experience and can therefore be regarded solely as guidelines. Our products are subject to continuous development. We therefore reserve the right to make possible changes.

General information:
Polymerization time for chemical polymerization approximately 3 min. (without light in the absence of air).
Optionally, the surfaces can be exposed to a suitable polymerization lamp for 20 sec.
Use a light polymerization system with an emission range of 350 - 500 nm to polymerize the material. The required physical properties are only reached if the polymerization light functions properly. Therefore, it is necessary to check the light intensity regularly according to the manufacturer's instructions.

Storage:
Store at 10 - 25°C (50 - 77°F).

Shelf-life:
The maximum shelf-life is printed on the label of each packaging. Do not use after the expiration date.

Side-effects:
With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us.

Contra-indications / interactions:
If a patient has known hyper-sensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. The dentist should consider known interactions and cross-reactions of the product with other materials already in the patient's mouth before using the product.

Note:
Please supply the attending physician with the above information, if this medical device is used to produce a special model. Please also consider the safety data sheets.

Disposal:
Must be specially treated adhering to official regulations.

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