**IMPLA**

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## **Implantat-System**

Artikel-Nr.: 608024

**SCHÜTZ DENTAL**  
Micerium Group

Schütz Dental GmbH • Dieselstr. 5-6 • 61191 Rosbach/Germany  
Tel.: + 49 (0) 6003 814-0 • Fax: + 49 (0) 6003 814-906  
[www.schuetz-dental.de](http://www.schuetz-dental.de) • [info@schuetz-dental.de](mailto:info@schuetz-dental.de)

**CE 0297**

Mandler 11/2023 – 500

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If there is anything in this instruction for use that you do not understand, please contact our customer service department before using the product.

As the manufacturer of this medical device, we inform our users and patients that all serious events occurring in connection with it must be reported to us (the manufacturers) as well as the relevant authorities in the Member State where the user and/or patient is resident.

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


## List of contents

Page

DIALOG Implant .....	3–7
General notes on implantation .....	8–9
IMPLA locking screws (covering and healing caps) and gingiva formers .....	10–11
IMPLA Abutments .....	12–13
Locator® Abutment – Abutment System .....	14–17
Multi-Unit Abutments .....	18–19
IMPLA PreFace Abutments .....	20
IMPLA Aesthura Abutments .....	21
Disinfection, cleaning and sterilization of surgical and prosthetic components of the IMPLA system .....	22–25
Application and safety instructions .....	26
Manufacturer’s information on rotating surgical instruments Standardized manual reprocessing of rotating instruments .....	27–28

# DIALOG Implant

IMPLA implants are made of pure titanium Grade 4 (ISO 5832-2). The IMPLA system comprises surgical, prosthetic and laboratory components and instruments.

## Important notes, warning

Implants are part of a complete system and may only be used with the original components and instruments of the system. The use of components which are not part of this system voids the guarantee and rules out parts being replaced. The IMPLA system may only be used by dentists, doctors and dental technicians familiar with dental implantology including diagnosis, preoperative planning, surgical techniques and prosthodontics.

Consultation is done in written form or by electronic media as well as with demonstrations during special workshops. All information is based on the state of the art and latest scientific findings known at the time of launching the product. However, we urge the users of our products to check the suitability of the product in each case. This is especially important for techniques not expressly recommended by the manufacturer. In case of doubt, we encourage our customers to confer with Schütz Dental Group/Section IMPLA.

## Reuse and sterilization

All abutments and components are for single use only and must not be reprocessed. Any reuse involves the risk of cross-contamination! All instruments, abutments and components are packaged non-sterile and must be sterilized (EN 13060) prior to placement (see pages 22 – 25). Using unsterile components increases the infection risk for the patient.

## Indications

DIALOG implants are available in different implant diameters and implant lengths:

Ø 3.7 mm	L: 8 / 9.5 / 11 / 13 mm
Ø 4.2 mm	L: 6.5 / 8 / 9.5 / 11 / 13 mm
Ø 4.8 mm	L: 6.5 / 8 / 9.5 / 11 / 13 mm
Ø 5.3 mm	L: 6.5 / 8 / 9.5 / 11 / 13 mm

DIALOG implants are used for:

- treatment of individual tooth gaps
- bridgework spanning interdental and free-end gaps according to general and special principles of prosthetic treatment guidelines
- restoration of edentulous mandibles and maxillae

### Further indications:

Implants with a diameter of 3.7 mm

These implants are an alternative in case of limited jaw ridge width. They are, because of their lower mechanical strength when compared with the standard implant (4.2 & 4.8 mm), intended for the following indications:

- a) replacement of individual mandibular incisors
- b) restoring lateral upper incisors
- c) four interforaminal implants to restore the edentulous mandible with a bar construction (primary splinting)

### Contraindications for implants with a diameter of 3.7 mm:

Implants with a diameter of 3.7 mm must not be used with ball anchors, magnets and telescopes.

Implants with a diameter of 4.2 mm

The standard implants with a diameter of 4.2 mm are indicated for all indication groups of oral endosseous implants in the maxilla and mandible. The prerequisite is an adequate minimum bone volume with a width of approx. 6.2 mm.

Implants with a diameter of 4.8 mm

The standard implants with a diameter of 4.8 mm are indicated for all indication groups of oral endosseous implants in the maxilla and mandible. The prerequisite is an adequate minimum bone volume with a width of approx. 6.8 mm.

Implants with a diameter of 5.3 mm

These implants are intended for use in the maxillary and mandibular posterior region (mainly molar region). The required alveolar ridge width should not be less than 7.3 mm.

The implants are primarily indicated for delayed implantations.

The implants with a length of 6.5 mm are generally not intended for single-tooth replacement or for accommodating a self-sufficient, prosthetic restoration. Implants with a length of 6.5 mm may only be used as a support implant, in conjunction (primary splinting) with implants > 6.5 mm in length.

### **Contraindications**

When selecting patients, the general contraindications for dental surgery must be taken into account. These include: Vascular diseases, uncontrolled diabetes mellitus and other metabolic or systemic diseases affecting bone regeneration, blood coagulation disorder, anticoagulative therapy, bone metabolism diseases.

Temporary contraindications: Acute and chronic inflammations, inadequate soft tissue present, chemotherapy or radiotherapy.

Oral contraindications:

- patient not really prepared to undergo total oral rehabilitation
- inadequate oral hygiene
- untreated parafunctions, e. g. bruxism
- inadequate amount of bone available
- inadequate interocclusal distance
- poor occlusion and/or articulation

### **Diagnostics / Patient information:**

Diagnostics consist of a medical history, clinical examination, X-ray examination using OPG and, if necessary, CT. A medical check up is recommended. It is also advisable to perform biomechanical treatment planning to avoid overloading the bone, the implant or the prosthetic construction, taking into account the individual masticatory forces and occlusion. Implantations are elective procedures that require extensive patient education. Treatment alternatives and consequences of the implantations must be pointed out. In particular, the importance of oral hygiene and the costs of follow-up treatment must be pointed out. With regard to the type and scope of the information, reference is made to the relevant case law.

### **Side effects:**

- Horizontal and vertical bone loss may come about even with a medically appropriate surgical procedure. The extent and the nature of bone loss cannot be estimated in advance.
- Reversible or irreversible impairments or interactions are possible if the implants are placed near special structures (e.g. adjacent teeth, maxillary sinus or nerves).
- In connection with dental implants, chronic pain may also occur in rare cases.

To minimize any risk of infection, antibiotic prophylaxis may be recommended. In rare cases, allergies to the prosthetic components (immediate implantation) can also prevent successful implantation.

### **Important notes**

When working with our implants and instruments intraorally, always use a protection against aspiration! Primary stability directly after the insertion of the implant is imperative for a successful healing. During the healing phase, IMPLA implants are generally not functionally loaded.

**Surgical procedure**

To create optimal conditions for the implant to heal, it is essential to treat the hard and soft tissues gently.

A special form drill is used for each respective implant configuration. This is matched to the implant geometry in terms of diameter and length. For use in hard and compact bone, the use of the cortical form drill is recommended to enable gentle, atraumatic insertion of the implant with the lowest possible insertion torque (maximum 30 Ncm). These serve to reduce the insertion torques for longer implant lengths. Since thermal trauma can hinder or prevent the healing of the implant, the temperature development should be kept as low as possible by the following measures:

- Recommended drill speeds for pilot drills/form drills/cortical drills:  
 $\varnothing$  1.8 mm: ca. 800 min<sup>-1</sup>,  $\varnothing$  2.2 mm: ca. 700 min<sup>-1</sup>,  $\varnothing$  2.7 mm: ca. 700 min<sup>-1</sup>,  $\varnothing$  3.0 mm: ca. 500 min<sup>-1</sup>,  
 $\varnothing$  3.7 mm,  $\varnothing$  4.2 mm,  $\varnothing$  4.8 mm,  $\varnothing$  5.3 mm: ca. 400 min<sup>-1</sup>.
- Intermittent drilling and milling technique with sharp drills and mills (max. 10 applications)
- Abundant external cooling (30 ml/minute) of the tool by physiological, sterile saline solution (NaCl) or Ringer's solution.

**Important**

To determine the optimum implant length and to prevent nerv injuries, check in advance that adequate space is available (X-ray).

**Insertion of the implant after preparation of the implant bed**

Depending on the design of the packaging, the implant can be inserted with the insertion aid and torque ratchet or with the implant key (insertion aid by machine) and contra-angle handpiece.

**Packaging with implant on the insertion post inside the glass tube:**

insertion post, IP, mount							
size	art./ref.	size	art./ref.	size	art./ref.	size	art./ref.
		4.2 / 6.5 mm	740185	4.8 / 6.5 mm	740191	5.3 / 6.5 mm	740197
3.7 / 8 mm	740180	4.2 / 8 mm	740186	4.8 / 8 mm	740192	5.3 / 8 mm	740198
3.7 / 9.5 mm	740181	4.2 / 9.5 mm	740187	4.8 / 9.5 mm	740193	5.3 / 9.5 mm	740199
3.7 / 11 mm	740182	4.2 / 11 mm	740188	4.8 / 11 mm	740194	5.3 / 11 mm	740200
3.7 / 13 mm	740183	4.2 / 13 mm	740189	4.8 / 13 mm	740195	5.3 / 13 mm	740201

**Insertion with insertion aid and torque ratchet**

After opening the blister seal, the implant is removed from the glass tube. For this purpose, the plastic plug is removed from the glass tube and the attached implant is carefully placed in the prepared implant site and gently screwed in until the implant carrier plug can be removed from the insertion post of the implant (if necessary with the aid of the open-end wrench). The insertion can now be continued – while maintaining the maximum insertion torque - using insertion tools for the torque ratchet.

**Removing the insertion aid**

In order to remove the insertion aid (for implants 3.7 mm brown; 4.2 mm red; 4.8 mm violet; 5.3 mm blue), pull off the plastic implant carrier and loosen the vertical screw (blue) with the hexagonal screwdriver SW 1.2 mm and dispose of the vertical screw. Secure the implant against rotation with help of a counter wrench in cases of low primary stability and concomitant augmentation in the according alveole. Afterwards, the insertion aid can be lifted from the implant.

**Packaging with implants without insertion aid inside an acrylic container:**

insertion post, IP, mount							
size	art./ref.	size	art./ref.	size	art./ref.	size	art./ref.
		4.2 / 6.5 mm	740485	4.8 / 6.5 mm	740491	5.3 / 6.5 mm	740497
3.7 / 8 mm	740480	4.2 / 8 mm	740486	4.8 / 8 mm	740492	5.3 / 8 mm	740498
3.7 / 9.5 mm	740481	4.2 / 9.5 mm	740487	4.8 / 9.5 mm	740493	5.3 / 9.5 mm	740499
3.7 / 11 mm	740482	4.2 / 11 mm	740488	4.8 / 11 mm	740494	5.3 / 11 mm	740500
3.7 / 13 mm	740483	4.2 / 13 mm	740489	4.8 / 13 mm	740495	5.3 / 13 mm	740501

**Insertion with insertion aid and torque ratchet**

After opening the blister seal, the implant is removed from the blister inside the container. The lid is moved to the side with the thumb, the container is now open. The implant is then removed directly from the container with a suitable insertion aid, which is already locked in the contra-angle handpiece. The container can be slightly compressed to improve the locking of the insertion aid in the implant. Now the implant can be inserted slowly into the previously prepared implant site in the jawbone with the contra-angle handpiece – while observing the maximum insertion torque.

**Soft tissue treatment**

Before wound closure, the cover and healing screw is screwed onto the implant. This is located at the top of the sealing plug of the glass tube in the case of packaging with insertion post or on the underside of the container in the case of packaging in a plastic container. The cover and healing cap indicate the prosthetic color code brown, red or blue.

Subsequent wound closure is performed with atraumatic suture material.

**Recovering the implant**

After a healing period to be individually determined, which should be no less than three months (mandible) or no less than six months (maxilla), the implant is exposed and the cover and healing screw are removed. If augmentative procedures are used, the healing period should be extended accordingly. The interior of the implant is then cleaned and a healing abutment is screwed in according to the thickness of the mucosa. Make sure that the gingiva former is seated correctly! A prosthetic restoration of the implant should only be performed after the soft tissue has healed without irritation.

**Taking the impression****Open technique**

After removing the gingiva former, a forming aid with an adequate diameter is attached to the implant with help of a fixing screw. If an impression is subsequently taken with a customized impression tray, the impression aid is first coated with impression material. After the impression material has set and the fixing screw has been loosened, the impression is removed together with the impression aid and the implant is closed again with the gingiva former. A laboratory implant is screwed to the impression post in advance and then a model is created – if possible with a mucosa mask. After the plaster has hardened, the impression aid or transfer abutment is removed and replaced by an appropriate prosthetic abutment.

**Closed technique**

Impression posts of suitable diameters are screwed onto the implant with help of a lab screw. Impressions are taken with a standard impression tray using transfer caps that are clicked onto the impression posts. The transfer caps remain within the impression and, after resetting the impression posts, are fitted with a laboratory implant to produce a plaster model.

**Laboratory components**

When selecting the abutments, the angulation of the implant and the thickness of the gingiva must be taken into account. Those areas of the abutment which are in contact with the implant may not be blasted with abrasives or treated in any other manner.

**Documentation**

Each IMPLA implant can be traced via its lot number which is noted on a self-adhesive label for internal user documentation. Possible claims can only be asserted by providing this data.

### **Packaging and sterility**

DIALOG implants are packed sterilely in a special clear container.

The selected implant is removed from the sterile container with the insertion aid only directly before insertion. It is inserted into the prepared bone cavity without direct contact with the implant surface using the corresponding insertion tool in the ratchet or contra-angle handpiece.

In case of packagings with an insertion post inside the glass cylinder, the covering and healing screw is screwed into the backside of the implant carrier. In case of packagings in an acrylic container, it is screwed into the bottom of the container. It can easily be retrieved with the hex screwdriver SW 1.2 mm.

Important reference data on implant length, implant diameter, article number, batch, CE certification, etc. are printed on the outer packaging made of cardboard.

### **General information**

- All implants are for single use only.
- Each is packed sterilely in a closed container, please note the expiration date of sterility.
- The implant is sterile (radio sterilisation) as long as the packing remains unopened and undamaged. A subsequent sterilization by the user is not permitted by the manufacturer.
- A chemical sterilization of the implant (chemiclave) must not be carried out under any circumstances, as this could damage the implant surface. In case of damaged packaging, the implant is no longer sterile and must not be used anymore.
- Do not touch the implant.
- Store the implant in a dry place at room temperature. Avoid direct light incidence.
- System errors due to mixing up instruments and implants should be avoided. Please always check the color codings and labelings.
- Thermal damages of the jaw bone by rotating instruments (such as pilot drills, thread cutters or enlargement drills) must be avoided.
- It is advisable that the implantation is performed only after the patient information consultation has been completed, after the patient has signed the consent form and after the patient has undergone a preliminary medical examination for implantation.  
We advise the implantologist to prepare a surgery protocol. It goes without saying that is necessary to perform implantological follow-up care in the sense of recall appointments to support the healing process and the patient behavior after implantation.

### **Torque of screw connections**

15 Ncm

- bar screw cone adapter
- Multi Unit secondary prosthetic screw

20 Ncm

- vertical screw
- cone adapter
- Multi Unit abutment 0°
- Multi Unit vertical screw
- Aesthura fixing screw

30 Ncm

- Locator abutment

## General notes on implantation

### Attention

All rotating instruments are subject to wear. Therefore, we cannot warrant the instruments' (pilot drills, extension drills, countersinks, cortical drills and thread cutters) functionality for more than 10 times of use, depending on the bone density. If the instruments are used more often, the manufacturer cannot guarantee a proper cutting efficiency and minimal temperature development any longer.

- Insufficient primary stability can be caused by poor preparation of the implant site or by insufficient bone quality. In this case, the loose implant should be removed and be replaced by an implant with a larger diameter. Implants should always be screwed into the bone all the way, even with their thread-free neck, when a covered healing process is considered (see above). Above all, a primary stability of the implant must be achieved.
- Implants should only be implanted when all prerequisites for a complication-free treatment and a longterm success are given. If there is any concern regarding the success of the implantation and/or the stability of the desired implant-based denture, a different solution without implantation should be favored.
- All rotating instruments can, in principle, be used several times. The user physician/expert alone decides under his/her own responsibility on the possibility of using the instrument again and the frequency of use on the basis of the respective case of use and the possible wear of the instruments. In case of doubt, the instruments should always be sorted out and replaced at an early stage.

### Care of instruments and parts

- Used instruments/parts must always be disinfected before cleaning.
- Damaged instruments/parts must be disinfected, cleaned, and disposed of.

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**Please consider our information on each material on page 22.**

### Disinfection

- Instruments/parts must be submerged in an adequate disinfection solution.
- Disinfection may be carried out manually or automatically. Please observe the manufacturer's instructions concerning dosing/concentration, residence time, temperature.
- In order to avoid damages, the instruments/parts may not have any contact during disinfection.
- It is highly recommended to wear protective gloves when working with contaminated instruments/parts (components).

### Cleaning

- Only use adequate cleaning agents.
- Always observe the manufacturer's instructions concerning concentration of the agent, temperature, etc.
- It is highly recommended to wear protective gloves (DIN EN 455) when cleaning contaminated instruments/parts manually.
- In case of automatic cleaning and disinfection, the instruments/parts may not have any contact in order to avoid damages.
- Only instruments/parts made of the same material should be cleaned together.

- In case of heavy contamination of the instruments/parts, ultrasound cleaning techniques are recommended.
- After cleaning, rinse the instruments/parts thoroughly with water and dry them immediately.
- Damaged instruments/parts must be removed and disposed of separately.

#### **Sterilization**

- All surgical instruments and parts must be sterilized after cleaning.

**Please consider our information on each material on page 22.**

#### **Storage**

- The sterilized instruments/parts must be stored in a dry and dust-free atmosphere at room temperature.
- After the expiration of the maximum storage period, the instruments/parts must be sterilized again.

Dose 25 kGy/min

# IMPLA locking screws (covering and healing caps) and gingiva formers

## 1. Indication

Surgical locking screws are used to immediately seal IMPLA implants after insertion and during submerged healing. The IMPLA gingiva formers can also be inserted during non-submerged healing or after successful submerged healing in order to form the gingiva surrounding the implant after the implant has healed.

## 2. Contraindications

None known

## 3. Notes

It is essential to read these instructions for use before using the IMPLA gingiva formers and locking screws.

When using the IMPLA gingiva formers and locking screws, the general rules for dental/surgical behaviour and indications, as well as the regulations for health and safety at work and the prevention of accidents must be taken into account. IMPLA gingiva formers and locking screws are only intended for medical/dental use in combination with the IMPLA implant system.

In the event of uncertainties in relation to the indication or type of application, use of the products should be ceased until all issues have been resolved. These instructions for IMPLA locking screws and gingivaformers are not sufficient to ensure proper use by handlers not experienced in implant-based prosthetics.

The user must be familiar with dental surgery and prosthetics, including diagnosis and preoperative planning and/or laboratory procedures. Since the use of this product is out of the control of the manufacturer, the latter shall not be liable in any way for damage caused. The responsibility rests exclusively with the user, or the attending physician.

When using IMPLA gingiva formers and locking screws, the following possible complications have been observed:

- The aspiration or swallowing of parts used in the patient's mouth
- Release or breakage of the screw thread due to incorrect torque (see recommended torque settings in Section 6 of the IMPLA instructions for use).

## 4. Precautionary measures

The following precautionary measures must always be taken before or during treatment.

- Before any procedures, ensure that all necessary parts, tools, and instruments are complete, functional, and available in the required amount.
- For your own safety, always wear suitable protective clothing.
- Position the patient so that the risk of aspirating the parts is minimised or excluded. Protect all parts used in the patient's mouth against aspiration and swallowing.
- Adhere to the specified torque for fixing prosthetic components (see recommended torque settings in Section 6 of the IMPLA instructions for use).

## 5. Side effects / Interactions

In very rare individual cases, allergies or sensitivities to the materials used to manufacture the prosthetic components (grade 5 titanium, etc.) cannot be ruled out.

In case of occlusal or approximal contact, the use of different alloy types in the same oral cavity may lead to galvanic irritation.

## 6. Use / insertion of locking screw / gingiva former

After fully inserting the implant in its final vertical position and loosening the insertion post, fit the surgical locking screw or gingiva former on the screwdriver (SW 1.2 mm), and screw into the implant by hand. The locking screws or gingiva formers are to be hand tightened (5–10 Ncm), taking into

account the primary stability of the implant. The gingiva former remains in for approx. 14 days to shape gingiva in place. The size of the gingiva former shall be selected according to the thickness of the gingiva, and the diameter and shape selected according to the prosthetic requirements.

**Note:** Gradually stretch to the desired diameter for very tight gingiva. Change to larger diameter every 5–7 days. Excessive pressure on the gingiva is to be avoided.

### **7. Cleaning / Sterilisation**

The IMPLA locking screws (included in the implant packaging) and IMPLA gingiva formers are provided sterile in a suitable blister packaging.

### **8. Warning**

The integrity of the product and sterile blister packaging must be checked before use. The product must not be used if the blister packaging is damaged.

## IMPLA Abutments\*

\*The term „abutments“ is used here to represent ball connector, adhesive base, cone connector, SAE cone adapter, IMPLA Direct base.

### **Important notes / warnings:**

IMPLA Abutments are a part of the IMPLA system and must only be used with the proper original IMPLA components. For the insertion of the abutments, only use the IMPLA instruments and tools. Otherwise, a damage of the original components cannot be ruled out.

### **Delivery:**

All IMPLA Abutments are supplied in unsterile packaging. Therefore, the components must be sterilized accordingly before use.

### **Material:**

IMPLA Abutments are manufactured from titanium grade 5 (Ti6Al4V ASTM 136).

IMPLA Scan abutments are made from PEEK (polyether ether keton).

### **Side-effects:**

Allergies and sensitivities towards the alloy cannot be ruled out completely and may occur in very rare cases.

### **Interaction:**

In case of occlusal contact, the presence of different alloy types in the patient's mouth may lead to galvanic reactions.

### **Safety advice:**

Metal dust and ceramic dust are hazardous to the health. When working on these materials, always wear a protective mask and protective goggles. Use an extraction system with a fine particle air filter commonly used in dental practices.

### **Indications:**

IMPLA ball abutments are used to anchor removable restorations to IMPLA implants. For this purpose, the prosthesis is ground out basally at a suitable area and a ball head matrix with an internal silicone O-ring is incorporated. The O-ring is subject to wear. It can be removed from the matrix and replaced, also for cleaning.

The IMPLA adhesive base is used to fix CAD constructed and CAM produced restorations on IMPLA implants. This requires a suitable design software as well as material blanks and a milling machine to produce the prosthetic structures, which can then be bonded to the adhesive base.

Cone connectors are intended for the production of cementable individual crowns and bridges.

SAE cone adapters are used for the production of bar abutments and bridges which are screw-retained in a tension-free method.

The IMPLA Direct base is intended for the construction of bar abutments or screw-retained bridges using the sprue technique, casting technique or laser technique.

Scan abutments are only intended for temporary use in the patient's mouth for the time of optical scanning using the intraoral scanner.

### **Contra-indications:**

When selecting the patients, it is extremely important to ensure that they are able to maintain the necessary oral hygiene. In addition, the general rules and limitations in surgery and prosthetics for the selection and treatment of patients apply.

### **Procedure:**

If the abutment does not bear a screw thread, carry out the final fixation on the intended implant with the vertical screw standard, art. no. 636649, unless otherwise specified for the respective abutment type in separate instructions.

For the following cone connectors, the final fixation on the intended implant may only be performed with the vertical screw standard, article 636649 from batch 2020000804.

- 636198 Conical connector 3.3 mm / 15° – GH = 3 mm
- 636183 Conical connector 4.2 mm / 15° – GH = 3 mm
- 636184 Conical connector 5.3 mm / 15° – GH = 3 mm
- 636186 Conical connector 4.2 mm / 20° – GH = 3 mm
- 636187 Conical connector 5.3 mm / 20° – GH = 3 mm
- 636190 Conical connector 3.3 mm / 0° – GH = 1 mm
- 636191 Conical connector 4.2 mm / 0° – GH = 1 mm
- 636192 Conical connector 5.3 mm / 0° – GH = 1 mm
- 636193 Conical connector 3.3 mm / 0° – GH = 3 mm
- 636194 Conical connector 4.2 mm / 0° – GH = 3 mm
- 636195 Conical connector 5.3 mm / 0° – GH = 3 mm

For the following cone connectors, the final fixation on the intended implant may only be performed with the vertical screw short, article 638634.

- 636185 Conical connector 3.3 mm / 15° – GH = 1 mm
- 636196 Conical connector 4.2 mm / 15° – GH = 1 mm
- 636197 Conical connector 5.3 mm / 15° – GH = 1 mm
- 636188 Conical connector 3.3 mm / 20° – GH = 1 mm
- 636189 Conical connector 4.2 mm / 20° – GH = 1 mm

## LOCATOR® ABUTMENT – Abutment System

**The LOCATOR® Implant Attachment System includes:** LOCATOR® Abutments, Retention Inserts, Denture Attachment Housing, Ancillary Processing Parts (including abutment analogs, processing spacer, block out spacer, parallelpost, impression coping, black processing male), and Tools (including LOCATOR® instrument, LOCATOR® adapter with connector for contra angle and corresponding adapter for torque ratchet).

**DESCRIPTION:** The LOCATOR® Implant Attachment System is a universal hinge, resilient attachment for endosseous implants in the mandible or maxilla in order to restore masticatory function. The attachment system allows for the prosthesis to be removed and replaced by the patient.

**INDICATIONS FOR USE** – The LOCATOR® Implant Attachment System is designed for use with overdentures or partial dentures, retained in whole or in part by endosseous implants in the mandible or maxilla.

**CONTRAINDICATIONS – LOCATOR® Abutment:** Not appropriate where a totally rigid connection is required. Use of a single implant with divergence of greater than 20 degrees from vertical is not recommended.

**MRI SAFETY INFORMATION** – A person with an implant or abutment device of the LOCATOR® Overdenture Implant System may be safely scanned under the following conditions. A scanning environment under Normal Operating Mode, with a Static Magnetic Field Strength ( $B_0$ ) of 1.5 T or 3.0T, with a Maximum Spatial Field Gradient of 40 T/m (4,000 gauss/cm), and RF Excitation which is Circularly Polarized (CP). The Maximum Whole-Body specific absorption rate (SAR) is 2 W/kg and the Maximum Head SAR is 3.2 W/kg, both under Normal Operating Mode. The Scan Duration limit is 2 W/kg whole-body average SAR for 60 minutes of continuous radio frequency (RF), (a sequence or back-to-back series/scan without breaks). Failure to follow these conditions may result in injury. NOTE: The presence of this device may produce an image artifact that can extend as far as approximately 23.0 mm  $\pm$  0.5 mm from the device, as observed in a nonclinical setting. The  $\pm$  0.5 mm is the numerical value of the combined standard uncertainty.

**STORAGE AND HANDLING** – The LOCATOR® Implant Attachment System in its undamaged, original packaging is not subject to any special considerations for storage or handling (during transport and storage).

**WARNINGS AND PRECAUTIONS** – The abutments of LOCATOR® Implant Attachment System are made from the titanium alloy Ti-6Al-4V and have a Titanium Nitride (TiN) coating. Restorative components are made from Ti-6Al-4V, 303 Stainless Steel and Nylon. Patients should be evaluated for allergies or hypersensitivity to such materials prior to treatment. During the treatment process patients may be exposed to instruments, tools, and processing components which may contain various grades of stainless steel, polyethylene, nitrile, or silicone.

Product should be inspected for integrity prior to use. Product from damaged packaging should not be used on patients. In the event that the packaging is damaged, the damaged packaging with the product should be returned to the manufacturer and a replacement will be provided only if damage to packaging is caused by product shipment.

If the LOCATOR® Implant Abutment is subjected to inappropriate loading conditions, there may be a potential risk of metal fatigue.

As surgical instruments are susceptible to damage and wear, they should be inspected before each use. Markings should be visible and legible. Any reusable instrument should be replaced if damage or wear is present to ensure proper functionality. The number of uses will vary and depends on a variety of factors including but not limited to bone density encountered, handling, proper cleaning, autoclave exposure, and storage conditions (do not store tools or instruments wet). Over time, repeat sterilization may affect appearance and visibility of markings. When applicable to the surgical instrument, check the connection feature for wear to ensure the connection is not damaged.

Patient evaluation including the determination of the general health, oral hygiene habits and status, motivation towards good dental care, and anatomical acceptability prior to the placement of the implant attachments as part of restorative process, is critical. Thorough evaluation of the patient's medical

status and health history is mandatory. Treatment planning is vital to the success of the implant and prosthesis.

The use of each of these attachment systems require that the clinician be thoroughly familiar with the product and the method for its use and application. The clinician must also use reasonable judgment in deciding when and where to use the product.

**SINGLE-USE DEVICES** – The LOCATOR® Implant Attachment System components, with the exception of the tools and instruments, are single-use devices and are provided non-sterile. Single-use devices must not be reused or re-sterilized. Reuse of a single-use device may cause harm to the patient in the transfer of blood, tissue, or saliva that may contain infectious disease. Single use devices may not function as intended if re-sterilized and may result in an improper surgical procedure and lead to improper function or failure of the device.

LOCATOR® Males: The inadvertent re-use of LOCATOR® nylon males could cause loss of retention of the overdenture due to wear from previous use or damage during removal with the LOCATOR® Core Tool.

LOCATOR® Attachments: The inadvertent re-use of LOCATOR® Attachments could contain patient contamination build-up and subsequent wear of the retention inserts. This would result in improper fit and function which cause the loss of retention of the prosthesis.

**MULTI-USE DEVICES** – The surgical instruments and tools of the LOCATOR® Implant Attachment System are multi-use devices. Reusable tools and instruments must be cleaned and sterilized prior to reuse on patients

TOOLS: The LOCATOR® Tools are designed for multiple uses, and are provided NON-STERILE. Follow the instructions provided here within for proper sterilization of non-sterile components, and the instructions for cleaning and reesterilization process of reusable components.

**CLEANING AND STERILIZATION** – The LOCATOR® Implant Attachment System, and other restorative components, instruments, and replacement LOCATOR® Attachments (sold separately), tools and instruments are supplied NON-STERILE, and should be sterilized prior to use on patients.

**CLEANING** – Reusable tools and instruments should be cleaned according to applicable instructions from the device manufacturer. **(i)** Disconnect or disassemble the instruments. **(ii)** Soak the instruments in enzymatic cleaning solution mixed according to manufacturer's instructions by completely submerging them for 20 minutes. Scrub the components using a soft-bristled, nylon brush until all soil has been removed. **(iii)** Remove the instruments from the enzymatic cleaning solution and rinse in tap water for a minimum of 3 minutes. Make sure to thoroughly flush internal holes/crevices of the instruments that have difficult to reach areas. **(iv)** Visually inspect instruments and tools for cleanliness and presence of residual debris. If additional cleaning is needed, place instruments in ultrasonic cleaner with enzymatic cleaning solution prepared according to manufacturer's instructions making sure that they are completely submerged, and sonicate for 10 minutes. **(v)** Remove the instruments from the ultrasonic cleaner, and rinse for 3 minutes making sure to thoroughly flush cleaning solution out of the holes/crevices and/or difficult to reach areas. **(vi)** Remove excess moisture from the instruments with a clean, absorbent, non-shedding wipe.

**STERILIZATION** – Tools and instruments provided non-sterile should be sterilized prior to use on patients. Reusable tools and instruments should be cleaned and sterilized prior to reuse on patients. The nylon males may be sterilized/disinfected using a liquid chemical sterilant as described below. Tools and instruments should be sterilized according to applicable instructions from the device manufacturer.

In order to ensure that the inserts are sterilized/disinfected, all microorganisms including Clostridium sporogenes and Bacillus subtilis spores are eliminated, the inserts must be soaked for a minimum of 3 hours in the liquid sterilant at room temperature. **Note:** An FDA approved liquid chemical sterilant for critical devices that are heat-sensitive and incompatible with sterilization methods such as steam and gas/vapor/plasma low temperature processes may be used following the manufacturer's directions for the sterilization of the device. Titanium abutments and stainless steel or other metal instruments may be sterilized by autoclave sterilization using the following parameters. For gravity cycle, place components in autoclave bag; and for pre-vacuum cycle, wrap the component with autoclave wrap

material and secure wrap with autoclave tape. Wrap the components using a wrap that is FDA-cleared for the indicated cycles.

Autoclave Sterilization Parameters are listed below:

Cycle Type	Description	Temperature	Exposure Time	Drying Time
Gravity	LOCATOR® Abutments, Tools and Instruments	132 °C / 270 °F	15 Min.	30 Min.
Pre-Vacuum		132 °C / 270 °F	4 Min.	20 Min.

Re-sterilizable instruments should be dried completely and stored in a clean and dry location at normal room temperature. Prior to instrument use, the exterior of any sterilized packaging should be inspected for integrity. Care must be exercised in the handling of wrapped or autoclave bagged instrument kits or instruments to prevent damage to the sterile barrier. If damage to the sterile barrier is observed, resterilization is recommended for reusable devices only. Single use devices should not be resterilized.

**DISPOSAL** – Dispose of used devices which pose a risk of infection according to facility clinical waste procedures and applicable local and state regulations.

**PROSTHETIC PROCEDURES** – Based on the results of the patient's pre-surgical assessment, the clinician should select and order the appropriate LOCATOR® Abutment based on the type of implant diameter, and tissue height.

It is imperative that all bone and soft tissue be removed from the superior aspect of the implant body to guarantee complete seating of the abutment.

Using a torque wrench, tighten the LOCATOR® Abutment to 30 Ncm. Use of higher torque values than recommended above could cause a fracture of the LOCATOR® Abutment.

**Impression and Stone Model Fabrication**

(i) With the LOCATOR® Abutments torqued in place, snap the impression copings on the abutments until they are seated firmly. (ii) Proceed by taking an impression. (iii) Remove the tray and snap an analog into each intaglio of the impression coping. (iv) Capture the abutment position in stone using standard methods for fabricating a laboratory stone model.

**Prosthesis Fabrication**

(i) Seat the LOCATOR® Retention Housings with the processing inserts on each of the abutments. (ii) Fabricate the prosthesis using standard laboratory techniques. (iii) When delivering the prosthesis, use the lowest retentive level insert to begin with and increase the retention if needed.

**Denture Attachment Housing Pickup Technique (Optional)** – (i) Place a block-out spacer around each abutment and press down.

(ii) Seat the LOCATOR® retention housings with the processing inserts on each of the abutments. (iii) Secure the retention housings to the prosthesis using auto-polymerizing or light cure acrylic or composite resin pickup technique.

**Prosthesis Delivery**

(i) Once the fit of the prosthesis is verified, remove the processing inserts from each retention housing using the core tool. (ii) Replace them with the lowest retention level inserts initially and increase the retention if needed. Firmly snap the prosthesis into place, ensuring that each insert is fully engaged onto each abutment.

**HEALING PHASE** – For delayed loading protocols: Relieve the denture to ensure the abutments are not in contact with any denture acrylic. A soft liner may be added to the denture to ensure patient comfort during the healing phase.

**PATIENT CARE** – Good oral hygiene is vital to attachment success. The patient should be made aware of the following:

(i) The LOCATOR® Abutments must be thoroughly cleaned each day to prevent plaque build-up and the patient should use a soft nylon bristle or end-tufted toothbrush with a non-abrasive toothpaste to clean the abutments. (ii) The coarse particles in abrasive toothpaste may scratch the surfaces of the abutments and cause plaque accumulation. (iii) An irrigation system is recommended to flush out debris from the inside of the LOCATOR® Inserts. (iv) The inserts are made of a soft plastic material (nylon) to allow the overdentures to be removed/replaced regularly. Plastic materials are subject to wear as part of normal use and may require replacement. (v) Bruxism wears the LOCATOR® attachments and may reduce the longevity of the retention inserts.

Patients should be instructed to maintain routine follow-up visits for hygiene and attachment function evaluation. Should a patient experience any discomfort or loss of retention of the overdenture, they should consult a dental professional.

Follow-up visits are recommended at 6 month intervals. Abutments must be re-tightened at follow-up visits to the torque specified above. Failure to re-tighten abutments could lead to screw loosening and abutment fracture. Patients should be examined for signs of inflammation around the implant abutments and for implant mobility.

**Inserting and Removing the overdentures** – The patient should be instructed on how to properly insert the overdenture. The patient should make sure they can feel that it is positioned over the abutments prior to applying pressure. The patient should use both hands and press down one each side and firmly snap the overdenture into place.

**NOTE:** The patient must not bite their overdentures into place as this force will result in improper wear of the abutments, including the retention inserts in the overdenture. Remove the overdenture by placing the thumbs under the edges of the overdenture flanges and pulling each side upward/downward simultaneously. Use of the tongue may aid in removal. Once removed, a thorough cleaning is recommended.


**Cleaning your Implant Retained Denture** – Instruct the patient to follow the protocol below to ensure the longevity of their overdenture. (i) Fill a washing basin with some warm water to prevent fracture of the overdenture. Apply non-abrasive toothpaste onto the soft bristle toothbrush and thoroughly clean every surface of the overdenture. (ii) Each night, remove overdenture and immerse in a cup of plain water.

**Further Information** – Traditional restorative protocols should be followed to process the attachments into the patient's overdenture. Standard overdenture care and maintenance should be followed in order to ensure the longevity of the restoration.

## Multi-Unit Abutments

### Description:

The IMPLA Multi-Unit Abutment Concept facilitates improved usage of the available bone structure. By means of an angulated insertion of the terminal implants, the available bone is used in distal direction which improves the stability. With the according angulation, the Multi-Unit Abutments ensure a uniform insertion direction of the prosthetics to facilitate a tension-free screw fixation. For this purpose, the IMPLA Multi-Unit Abutment Concept contains straight and angulated components.

### Components Multi-Unit System

Abutments with the diameters 3.3 mm, 4.2 mm, 5.3 mm, each with gingiva heights of 1 mm and 3 mm, and angulations of 0°, 20° and 30° (not 5.3 mm).

### Accessories IMPLA Multi-Unit System

Impression system with impression aid and fixing screw as well as the according lab implant.

Gingiva sleeve for conditioning of the mucosa as well as the appropriate metal or plastic sleeve for the prosthetics with suitable screw for secondary screw fixation.

### Important notes / warnings:

IMPLA Multi-Unit Abutments are a part of the IMPLA system and must only be used with the proper original IMPLA components (IMPLA implants). For the insertion of the abutments, use only the IMPLA instruments and tools. Otherwise, a damage of the original components cannot be ruled out.

### Indications:

The IMPLA Multi-Unit Abutments are exclusively meant for a use with occlusal screw-fixed prosthetics (e. g. bridges, hybrid prostheses and cast bars).

A minimal implant length of 11.5 mm in the distal position is indicated for the use of IMPLA Multi-Unit Abutments.

The minimal number of implants for a complete reconstruction is four in the mandible and six in the maxilla.

A sufficient bone quality is a basic requirement for implantation, as it ensures the primary stability. The proper torque is between 30 and 40 Ncm.

The distal implants must be inserted with an angulation that corresponds to the angulation of the abutments inserted later, i.e. **0°, 20° or 30°**.

Presurgical planning as well as a comprehensive medical history and diagnosis (including radiological examination) are absolutely necessary.

In principle, the general rules and instructions for surgical and prosthetical procedures in dental implantology apply.

**Contra-indications:**

When selecting the patients, it is extremely important to ensure that they are able to maintain the necessary oral hygiene. In addition, the general rules and limitations in surgery and prosthetics for the selection and treatment of patients apply.

**Delivery:**

All components of the IMPLA Multi-Unit Abutment Concept are supplied in unsterile packaging. Therefore, the abutments must be sterilized accordingly before use.

**Procedure:**

Based upon presurgical planning and following the drill protocol for IMPLA implants, a pilot drilling for a mesial implant without angulation is performed first.

Subsequently, all other implants can be aligned to this first drilling with help of a drilling template. The drilling template improves the spacial orientation and presets the angulation of the distal implants. The drill is guided through a sleeve in the drilling template into the bone in the exact angle intended.

After the cavities have been prepared in accordance with the IMPLA drilling protocol, the implants can be inserted.

The prosthetic restoration can be carried out either in two stages or by immediate load.

In case of a two-stage restoration, we recommend to cover the implant with a surgical cover screw (covered healing) or with a gingiva former for open healing.

A primary stability of 30–40 Ncm is absolutely necessary for the immediate load approach.

Impressions are taken with help of the according impression components of the IMPLA Multi-Unit Abutment Concept.

Depending on the previously selected implant angulation (0°, 20° or 30°), the mesially positioned implants are then restored with a straight abutment and the distally positioned implants with an angulated abutment. The abutments are fixed with the appropriate vertical screw with a torque of 20 Ncm.

Finally, the prosthetic construction is set upon the abutments tension-free and is fixed with a prosthetic screw with a torque of 15 Ncm.

**Material:**

The IMPLA Multi-Unit Abutments are made of titanium grade 5 (Ti6Al4V ASTM 136).

**Side effects:**

Allergies and sensitivities towards the alloy cannot be ruled out completely and may occur in very rare cases.

**Interaction:**

In case of occlusal contact, the presence of different alloy types in the patient's mouth may lead to galvanic reactions.

# IMPLA PreFace Abutments

## Description:

The IMPLA PreFace Abutment is finished with an appropriate milling system with help of the original PreFace Abutment Holder. The PreFace Abutment Holder serves to fix the PreFace Abutment inside the milling system in the appropriate working position. The PreFace Abutments are held in the PreFace Abutment Holder on the front side, making sure that they are seated exactly in the end position. The PreFace Abutment Holders are available for a wide range of milling machines.

## Important notes / warnings:

IMPLA PreFace Abutments are part of the IMPLA system and must only be used with original components of the system (IMPLA implants). Process the PreFace Abutments only with the corresponding suitable abutment holders.

## Indications:

The IMPLA PreFace Abutments are used for the production of individual titanium abutments on DIALOG implants in CAD/CAM milling machines. The individual abutments can be combined with crowns or bridges to reconstruct function and aesthetics.

## Contraindications:

The IMPLA PreFace Abutment must only be combined with DIALOG implants. Manipulating the connector geometry will lead to an imprecise fit which will render the abutment useless.

## Delivery:

All components of the IMPLA PreFace Abutment concept are supplied in unsterile packaging. Therefore, the components must be sterilized accordingly before use by the practice staff.

## Procedure:

The shape of the PreFace Abutment can be adjusted individually to the present anatomic conditions. The production of the PreFace Abutment is done digitally with an appropriate CAD software. The basic rules must be adhered to.

Design the emergence profile gently. Any further expansion of the gingiva than the expansion obtained with the gingiva former should always be discussed with the attending dentist beforehand. The gingiva can be expanded gradually by individually designed temporary abutments.

In general: The circular notch should be positioned slightly towards the subgingival direction in the vestibular area and orally on the gingiva line. This way, the cement between coping and abutment is easy to remove and the transition is not visible from vestibular direction. The minimal wall thickness of the abutments should not fall below 0.5 mm (screw hole to the outer surface of the abutment). The diameter of the emergence profile should not exceed 13 mm, while the height of the emergence profile should not exceed 6 mm. The total length of the PreFace Abutment should not exceed 15 mm, the maximum angle of the PreFace Abutment should not exceed 30°. The taper of the PreFace Abutment should not exceed 6°. The PreFace Abutment should be rounded occlusally and should taper off. Avoid sharp edges.

The recommended torque for the fixation of IMPLA PreFace Abutments is 20 Ncm.

## Finishing:

In case that further work steps are necessary after shaping the PreFace Abutment in the CAD/CAM milling machine, please note the following:

To mill PreFace Abutments, use appropriate, flawless mills with only little pressure. The minimal wall thickness must not fall below 0.5 mm. Avoid ridges and edges. We recommend to fix the PreFace Abutments in the lab analog to protect the connection geometry.

## Scan body:

A scan abutment is used to determine the exact position of the implant during the scanning procedure. The scan abutment's diameter must correspond to the diameter of the selected implant.

\* PreFace is a registered trademark of Medentika GmbH

## IMPLA Aesthura Abutments

**Description:**

The IMPLA Aesthura Abutment is used to fabricate occlusally screw-retained restorations. It is screwed onto the intended implants with the vertical screw, art. no. 638665.

**Important information / warning:**

The Aesthura Abutment must not be altered by measures of any kind.

IMPLA Aesthura Abutments are part of the IMPLA system and must only be used with original components of the system (IMPLA implants). For the insertion of the abutments, use only the IMPLA instruments and tools. Otherwise, a damage of the original components cannot be ruled out.

**Indications:**

Individual gaps, partially or completely edentulous jaws, screw-fixed individual tooth restorations.

**Contra-indications:**

When selecting the patients, it is extremely important to ensure that they are able to maintain the necessary oral hygiene. In addition, the general rules and limitations in surgery and prosthetics for the selection and treatment of patients apply.

**Delivery:**

All IMPLA Abutments are supplied in unsterile packaging. Therefore, the components must be sterilized accordingly before use by the practice staff.

**Material:**

The IMPLA Aesthura Abutments are made of grade 5 titanium (Ti6Al4V ASTM 136).

**Side effects:**

Allergies and sensitivities towards the alloy cannot be ruled out completely and may occur in very rare cases.

**Interaction:**

In case of occlusal contact, the presence of different alloy types in the patient's mouth may lead to galvanic reactions.

**Procedure:**

The abutments are fixed with the corresponding Aesthura vertical screw at 20 Ncm.

## Disinfection, cleaning and sterilization of surgical and prosthetic components of the IMPLA system

Following surgical and prosthetic instruments/tools of the IMPLA system are composed of surgical steel/titanium/PEEK/PP and are to be treated accordingly:

X-ray spheres, initial drill, taper drill, taper cortical drill, GS initial drill, GS taper drill, GS taper cortical drill, GS anchor drill, drill extensions, ratchet, screw drivers, depth gauge, sealing and healing cap, gingiva former, abutments, scan abutments, impression aid, vertical screw, bar screw, bone screws, DIALOG tray, DIALOG GS tray, GS adapter, GS anchor screw

<b>Please observe!</b>	Neutralize and dry the instruments thoroughly in order to avoid precipitation on the instruments to be sterilized. It is recommended to autoclave unmixed materials only.
Disinfection	Always wear protective gloves when disinfecting components suited to this treatment. Avoid drying of sticking residues of blood and/or tissue. Always use an appropriate disinfection bath and dose according to the manufacturer's instruction.
Cleaning	The components can be cleaned with all cleaning agents suitable for reprocessing surgical instruments and are resistant to any cleaning process, e.g. ultrasonic bath. Please avoid the use of metal brushes and abrasives. After cleaning, subsequently clean pH-neutrally and rinse off the cleaner and particle residues with demineralized water. Dry immediately.
Sterilization	The components are resistant to any kind of sterilization process. Make sure that all deposits have been removed by cleaning. The parts may not have any contact during sterilization. The temperature should not exceed 176°C (248°F). Please observe the instructions of the autoclave manufacturer.

**Instructions for the practice:**

1. Do not let surgical residues (blood, tissue, etc.) dry on the instruments/parts. Submerge the instruments/parts in disinfective agent directly after the surgery.
2. Clean the instruments/parts only with a nylon brush. (Do not use metal brushes, such as brass brushes, etc.)
3. Rinse off disinfectant residue thoroughly with clear water.
4. Dry instruments/parts directly after cleaning and rinsing.
5. Never sterilize instruments made of different materials (e.g. ratchet) in the same cycle!
6. Please observe the instructions for use of the cleaning agents, as suitability for all instruments is limited (see material groups).

**1 Use of the torque ratchet:**

The torque ratchet is build up of the ratchet head (Fig. 1-b), the ratchet handle (Fig. 1-c) and the setting screw (Fig. 1-e). Inside the ratchet head, there is a tool fitting (Fig. 1-a).


Fig. 1

The ratchet head is lockable to lock the torque limit (see 1.1 below). If the ratchet head is not locked, the torque can be adjusted infinitely between 10 and 30 Ncm with the setting screw (see 1.2).

Before its first use and after each use in a patient case, the ratchet must be disassembled, cleaned and dried. **It must be sterilized directly before each use** (see section 3).


Fig. 2


Fig. 3

Fig. 4


Fig. 5

**1.1 Adjusting the ratchet**

In the factory setting, the ratchet head is not locked. If you would like to lock it, turn the setting screw to less than 10 Ncm (Fig. 4), to relieve all pressure on the spring.

Subsequently, pull the handle from the ratchet head (Fig. 2). Then, turn the handle clockwise by 90° and let it lock into place in the ratchet head. Finally, preset a torque of 20 Ncm to ensure that the ratchet remains locked.

**Attention:** In the fixed condition, the maximum load is 100 Ncm!

**1.2 Setting the torque function**

First, ensure that the ratchet is not locked. Turn the setting screw until you have reached the necessary torque. The correct torque is set once the according marking ring at the end of the handle is just still visible (Fig. 4). The ratchet head will now bend once the preset torque has been reached (Fig. 5). Afterwards, do not keep on turning, as this could lead to exceeding the preset torque significantly.

**Important notice:**

Whenever the ratchet is not in use, turn back the setting screw to relieve all pressure from the spring.

**1.3 Change tool fitting**


(See 2.1 [1.] and 2.2 [5.]

## 2 Maintenance of the torque ratchet

### 2.1 Disassembly

1. The tool fitting (Fig. 6-a) is removed from the ratchet head (Fig. 6-b) by pulling on the two pin ends (arrows) towards the ratchet handle. Please note that the tool fitting can only be pulled from the ratchet head towards one side. **Do not use force.**
2. Completely screw off the setting screw (Fig. 6-e) from the ratchet handle (Fig. 6-c).
3. Pull off the ratchet handle (Fig. 6-c) from the ratchet head (Fig. 6-b) and remove the spring.
4. Disinfect, clean and dry all parts (see section 3).

Afterwards and **before the sterilization**, the ratchet is reassembled.

### 2.2 Reassembly

1. Before reassembling the ratchet, all parts must be disinfected, cleaned and dried completely. Cover the areas marked with ● with oil (see 3.4).
2. Push the ratchet handle (Fig. 6-c) over the threaded bolt into the ratchet head (Fig. 6-b).
3. Push the spring (Fig. 6-d) into the ratchet handle (Fig. 6-c).
4. Screw the setting screw into the ratchet handle (Fig. 6-e) until the resistance of the spring becomes evident.
5. The tool fitting (Fig. 6-a) is set back into the ratchet head (Fig. 6-b) by pulling the two pin ends (arrows) towards the ratchet handle. Please consider that the tool fitting can only be set into the ratchet head from one side. **Do not use force.**

#### Important information:

Always relieve pressure on the spring when ratchet is not in use. Screw off the setting screw until the pressure on the spring is relieved completely.

## 3 Reconditioning

### 3.1 Disinfection

Before first use and immediately after each use, all ratchet parts must be disinfected in a suitable disinfectant. Adhere to the instructions of the disinfectant's manufacturer.

Chemical as well as thermal disinfection techniques are suitable, if they come up to the requirements of medical instrument disinfection.

### 3.2 Cleaning

Immediately after disinfection, all parts of the ratchet must be cleaned with a suitable cleaning agent. Adhere to the instructions of the cleaning agent's manufacturer. Otherwise, impurities will dry off onto the instrument and make cleaning difficult.

In order to avoid corrosion and discoloration, closely adhere to the instructions on concentration, submergence time and temperature given by the manufacturer of the cleaning agent.

Remove any deposits, impurities, etc. thoroughly with a soft brush (do not use metal brushes, scouring agents, wire wool or too much pressure).

#### Please note:

When disinfecting or cleaning the device with an ultrasound unit, please adhere to the instructions of the unit's manufacturer.

### 3.3 Rinsing and drying

Directly after disinfection, cleaning or ultrasound treatment, the individual parts must be rinsed carefully with clear desalinated or distilled water (do not use tap water). Afterwards, dry all parts thoroughly to avoid water stains.

If the individual parts are not rinsed, residue is not removed and/or the parts are not dried completely, this could lead to discoloration or corrosion of the device during sterilization.

### 3.4 Care and checking

Parts that have joints must be treated with a sterilizable oil spray (paraffin-based) after all cleaning and disinfection procedures. Such a spray will reduce metal-on-metal-friction, avoids fretting corrosion and reduces the danger of oxidation during heavy use.

### 3.5 Sterilization

Only sterilize carefully cleaned, dried devices without any type of corrosion.

Badly cleaned and insufficiently dried parts bear the danger of burn-in or corrosion damage.

If corroded devices are sterilized with immaculate instruments, flash rust can cause damages on the immaculate instruments.

Please note the instructions of the manufacturer.

#### Note:

Only sterilize the torque ratchet directly before use.


### 3.6 Storage

#### Important information regarding the storage of the ratchet:

The spring of the ratchet should be fully relieved immediately after use to ensure proper use in the long term. Only store the ratchet in a relieved state.

## Cleaning and sterilization of the ratchet

Screw off turn-knob, remove spring and handle.

To clean, pull back the interlock to unlock the receptacle

Sterilization of the ratched only in disassembled state, without the spring.  
Attention: Autoclave spring separately!

## Application and safety instructions

All rotating instruments and parts for medical use (dental implantology, surgery) were developed for their specific application. For that reason, incorrect use may cause premature wear of the instruments and endanger the patients.

Rotating instruments for the medical field must only be used by physicians and other experts who are intensively familiar with the use of these instruments due to their training and experience.

### Correct application

- Take care to use only technically and hygienically flawless, well-maintained and cleaned contra-angle handpieces.
- The instruments must be inserted correctly.
- The instrument must revolve at operating speed before being placed on the object (see page 5).
- Avoid wedging or levering the instrument as this may increase the risk of fracture.
- Thermal damage through the use of rotating instruments must be avoided (we recommend the slow-speed mode with sufficient cooling). Please see Surgical Procedure on page 5 of this manual.

### Pressing forces

- Excessive pressing forces must be avoided at all times. Excessive pressing forces may damage the working sections of the instruments or fracture their blades. This also results in an increased generation of heat.

### Cooling

- To prevent excessive heat generation during preparation, provide for adequate air or water cooling (at least 50 ml/min).
- Inadequate cooling may lead to irreversible damages of the bone and/or tissue.

## Manufacturer's information on rotating surgical instruments

### Product:

The present manufacturer's information applies to all instruments that are used for surgical application. These are rotating instruments made of stainless steel.

**Instruments delivered in non-sterile condition must be prepared prior to first use.**

### Limited number of reprocessing cycles:

The end of product life is basically determined by wear and damage from use. If necessary, known restrictions on the frequency of use of instruments must be observed (see page 5, „Surgical procedure“). Frequent reprocessing does not affect the performance of these instruments.

### Work station:

Hygienic precautions according to the provisions valid in your country.

### Storage and transport:

Place instruments in a cleaning/disinfection tank filled with a suitable detergent/disinfectant immediately after use in the mouth to prevent drying of residues on the instruments (protein fixation). It is recommended to process the instruments within one hour of use at the very latest. The instruments should be transported to the reprocessing site in the disinfection tank. Endodontic instruments can be transported in a special interim support equipped with a foam insert drenched in disinfecting solution.

## Standardized manual reprocessing of rotating instruments

### Equipment used:

- Nylon brush
- Suitable detergent/disinfectant for rotary instruments with proven disinfecting effect
- Ultrasonic bath

### Reprocessing:

1. Rinse off surface contamination under running water. Remove stubborn contamination with a nylon brush under running water, turning the instrument constantly.
2. Rinse instrument thoroughly with running water.
3. Place the instruments in a suitable sieve into the ultrasonic device filled with detergent/disinfectant.
4. Observe the instructions of the manufacturer regarding concentration and immersion time for chemical cleaning/disinfection in the ultrasonic device. Be sure to observe the full correct immersion time which does not start until the last instrument has been placed into the ultrasonic device. Attention: Do not exceed 45°C (113°F) (risk of protein coagulation)!
5. On completion of the immersion time, rinse instruments thoroughly with suitable water (preferably with demineralised water to avoid residues).
6. Dry instruments (according to RKI recommendations, preferably with compressed air).
7. Visual examination with a suitable magnifying device to ensure that the instrument is clean and undamaged (experience has shown that a magnification factor of 8 permits a visual examination). If there are still visible residues of contamination, repeat the cleaning and chemical disinfecting process until no visible contamination is left.

**Control and functional test:**

Instruments showing the following defects are to be discarded immediately:

- Blunt and chipped blades
- Deformations (e. g. bent/twisted/fractured instruments)
- Corroded surfaces

**Packing:**

Make sure that the packaging is suitable for the instrument and the chosen method of sterilization. Single pack: The packaging must be large enough to ensure that there is no pressure on the seal. In the set: Place instruments onto the tray provided or onto universal sterilization trays. The instruments must be protected. Use an appropriate method to pack the tray. Instruments with limited times of use are to be marked accordingly.

**Sterilization:**

Steam sterilization using a vacuum process at 134°C (273°F) in a device that complies with the provisions of DIN EN 13060; with validated processes.

- Fractionated pre-vacuum (type B) respectively simplified pre-vacuum (type S).
- Sterilization temperature: 134°C (273°F)
- Hold time: 5 minutes minimum (full cycle)
- Drying time: 10 minutes minimum

In order to prevent staining and corrosion, the steam must be free of particles. The recommended limits for particle contents in feed water and condensed steam are defined by standard DIN EN 13060. Make sure not to exceed the maximum capacity of the sterilizer when sterilizing several instruments. Follow the instructions of the device manufacturer.

**Transport and storage:**

The packed sterile goods must be protected from dust, moisture and recontamination during transport and storage.

**General note:**

Observe the legal provisions regarding the reprocessing of medical products valid in your country (e. g. [https://www.rki.de/EN/Home/homepage\\_node.html](https://www.rki.de/EN/Home/homepage_node.html)).

The manufacturer confirms that the above detailed reprocessing methods are **suitable** for preparing the above named instrument group to enable their reuse. The person carrying out the reprocessing is responsible for ensuring that the applied method is carried out with appropriate equipment, materials and trained personnel at the reprocessing site and that it actually achieves the desired result. To guarantee this, routine controls of the validated mechanical and/or manual preparation methods are necessary. Any deviation from the above detailed process (e. g. use of different chemicals) must be carefully checked by the operator to ensure effectiveness and to avoid possible adverse consequences.