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If there is anything in this instruction for use that you do not understand, please contact our customer service department before using the product.

As the manufacturer of this medical device, we inform our users and patients that all serious events occurring in connection with it must be reported to us (the manufacturers) as well as the relevant authorities in the Member State where the user and/or patient is resident.

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#### Ceramic-Composite Repair System

##### Indication

The combination of the ReVeneer Primer and Base represents a fast and universally applicable system for intraoral and extraoral repair of fractured ceramic veneers with acid etching.

##### Preparation of Surface Fractures

Smooth the damaged surface and roughen it with a fine grit diamond (< 50 µm) or intraoral air polisher before cleaning off any adhering material with oil-free compressed air. If an appropriate moisture control cannot be achieved with cotton rolls, apply a rubber dam.

##### Please note

The use of an intraoral air polisher (Al<sub>2</sub>O<sub>3</sub> / < 120 µm) further enhances the bond strength.

ReVeneer Primer
<b>Application</b>
In combination with ReVeneer Base, ReVeneer Primer serves as an ahesive agent for ceramic /metal and composites.
Using a brush, apply a thin coat of primer to the ceramic or any exposed metal immediately after cleaning with oil-free compressed air.Allow to air-dry for approximately 1 minute. A second coat may be applied if necessary. The treated surface must be completely dry before continuing, otherwise adhesion will be limited significantly. Do not touch the surface after the primer has been applied.
<b>ReVeneer Opaque</b>
<b>Application</b>
ReVeneer Opaque serves as a coating for the metal frame. Its use is only necessary if the fracture extends to the metal frame.
Using a hard, shorthaired brush, apply a thin coat of the opaque. An insuffi-ciently opaque coat adversely influences the color effect. At least two coats are necessary to achieve total coverage. The fluidity of the opaque is opti-mized through stirring on a mixing block. Polymerize each layer separately with a conventional polymerization lamp (for polymerization times, please see table 4). For polymerization, hold the light guide as close to the surface as possible. Do not touch the surface of an opaquer layer.
Thick opaquer layers prevent an optimal polymerization and thus weaken the bond. An extension of the polymerization time does not result in a greater curing depth.
<b>Processing range of the mixture</b>
1–3 minutes, depending on light conditions.
<b>ReVeneer Base</b>
<b>Application</b>
ReVeneer Base is not intended for a completion of the veneer. It solely improves the adhesive bond between framework material and filling composite. The actual repair is carried out subsequently with a conventional micro hybrid composite. <b>Please note:</b> As many nano hybrid composites feature a high E-module, we discourage the use of such filling composites.
<b>Application</b>
Apply ReVeneer Base to the surface fracture with the supplied curved applica-tion tip. Spread the material evenly over the complete surface with a brush.
For hygiene reasons, the supplied brush inserts and application tips are disposable and are meant for one use only.
Polymerize the applied layer with a conventional polymerization lamp (for polymerization times, please see table 4). Hold the light guide as close to the surface as possible. During polymerization, a dispersion layer is formed on the surface which must not be touched or removed. It serves as a bond to the following layer.
<b>Polymerization times for ReVeneer Opaque and ReVeneer Base</b>
<small>(Errors and omissions excepted.)</small>
<b>Polymerization of ReVeneer Opaque and ReVeneer Base</b>
For the polymerization, use a light polymerization lamp with an emission range of 350–500 nm. The required physical properties are only obtained with properly functioning lamps. Regularly check the light intensity according to the manufacturer’s instructions.
<b>Intraoral Polymerization, Chair-side</b>
• Type of polymerization lamp: LED or halogen lamp
• Light intensity: ≥ 650 mW/cm²
• Wavelength: 350–500 nm
• Polymerization time ReVeneer Opaquer: 60 sec. per layer
• Polymerization time ReVeneer Base: 40 sec. per layer

Use conventional light polymerization lamps for the extraoral polymerization at the lab, e.g.:

Polymerization lamp	ReVeneer Opaquer	ReVeneer Base
Spektra LED	60 sec.	20 sec.
Spektra 2000 fast	180 sec.	180 sec.
Uni-XS / HiLite	90 sec.	90 sec.
GC Labolight LV-II / III	60 sec.	20 sec.
Solidilite EX / V	60 sec.	20 sec.

##### Please note

When creating large and time-consuming restorations either cover the material with a lightproof foil or temporarily move the operating light away from the site in order to avoid premature polymerization of the material.

##### Filling composite

Apply, polymerize and finish the composite according to the manufacturer’s instructions.

**Please note:** As many nano hybrid composites feature a high E-module, we discourage the use of such filling composites.

##### Important!

Please ensure that the treated surface is not contaminated by saliva or other contaminants during the complete repair procedure. Such contamination would limit the adhesive bond significantly. Avoid any contact between ReVeneer Primer and Base with the mucous membranes. If this condition cannot be provided with help of cotton rolls, a rubber dam must be applied.

This product was developed specifically for the described range of applications and must only be used accordingly. The manufacturer is not liable for damage caused by improper handling or processing.

##### Disinfection/ Protection from cross-contamination ReVeneer Base

Place the functioning syringe with attached application tip into a suitably shaped barrier sheath; pierce end of sheath with cannula, exposing the cannula for use. Using a barrier sheath facilitates cleaning and disinfection of the syringe between patients. After use of sheathed syringe, remove delivery tip and sheath by grasping on the hub of the delivery tip through the sheath, twist and remove tip along with sheath. Discard used tip and sheath in appropriate waste stream. Replace syringe storage cap. Disinfect – After removal of the delivery tip and sheath, disinfect this product using an inter-mediate level disinfection process (liquid contact) as recommended by the Centers for Disease Control and endorsed by the American Dental Asso-

ciation, Guidelines for Infection Control in Dental Health-Care Settings – 2003 (Vol. 52; No. RR-17), Centers for Disease Control and Prevention (USA).

##### Brush for Primer and Opaquer application

When using the supplied brush holder, please proceed as follows:

Place the brush holder into a suitably shaped barrier sheath; pierce end of sheath with the angled brush holder tip. Set the brush insert into the holder, exposing the brush for use. Using a barrier sheath facilitates cleaning and dis-infection of the brush holder between patients. After use of sheathed brush holder, remove brush insert and sheath by grasping the brush insert through the sheath, twist and remove it along with sheath. Discard used brush insert and sheath in appropriate waste stream. Disinfect – After removal of the brush insert and sheath, disinfect the brush holder using an intermediate level disinfection process (liquid contact) as recommended by the Centers for Disease Control and endorsed by the American Dental Association, Guidelines for Infection Control in Dental Health-Care Settings – 2003

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