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If there is anything in this instruction for use that you do not understand, please contact our customer service department before using the product. As the manufacturer of this medical device, we inform our users and patients that all serious events occurring in connection with it must be reported to us (the manufacturers) as well as the relevant authorities in the Member State where the user and/or patient is resident.

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## Instructions for use

# Sebond Grip

**Designation of purpose:** Framework covering materials are light-curing or chemically curing materials on composite basis for covering different framework materials for subsequent veneering with a veneering composite or a denture acrylic.

**Brief product description and users: Product description:** Light-curing universal groundig opaquer for the use between Sebond Smart and conventional light-curing powder or paste opaques which are available in the market.

Sebond Grip improves the chemical bond between primer and opaquer and contributes to an easy, safe and lasting bond.

**Users:** For use by dental technicians at a dental laboratory.

**Composition:** Silicon dioxide, urethane dimethacrylate, glass powder, pigments, initiators. Filler content: 62 weight %, anorganic filler (0.005 - 3.0 µm).

**Indications:** For covering framework materials in the crown and bridge technique (conventional dental metals or alloys) when producing new restorations and for repairs.

**Contra-indications:** Do not use in case of a known allergy against one of the components.

**Additional information:** Contains diphenyl(2,4,6-trimethylbenzoyl)-phosphine oxide. May cause allergic skin reactions.

**Cross-reactions with other agents:** – Phenolic substances (such as e. g. Eugenol) inhibit polymerization. Do not use materials which contain such substances. – Generally known cross-reactions or interactions of the medical device with other materials already in the mouth must be considered by the dentist when using the material.

**Framework design and preparation: Preparatory tasks:** The frameworks are designed, cast or CAD/CAM-milled and then finished as usual, in accordance with the standards of dental technology.

Conventional dental metals and alloys, PMMA, PEEK, and zirconium dioxide can be used as framework materials. All materials are to be used and to be prepared for processing (sand-blasting, conditioning) in accordance with the manufacturer's instructions. Before applying Sebond Grip, the frameworks must be conditioned. The conditioning of a framework must be carried out depending on the framework material used. Please observe the manufacturer's instructions. The surface is sand-blasted and blown off with oil-free compressed air.

As bonding agent (primer), use Sebond Smart in accordance with the manufacturer's instructions.

**Procedure:** Thinly apply Sebond Grip with help of a tight short-haired brush. At least two layers are necessary in order to achieve an optimum bond. It is possible to optimize the flowability by stirring on a mixing block. Make sure to apply an extra-thin layer in the area of retentions. Each layer is polymerized separately.

Light curing system	Time
HiLite / Uni-XS	90 sec.
Labolight LV-II/ III	1 min.
Solidilite	1 min.

**Subsequent tasks:** Afterwards, a conventional opaquer is applied and the restoration is designed with a conventional veneering composite. Observe the manufacturer's instructions.

**Repairs:** Carry out any repairs outside the mouth. The surface to be repaired is prepared as described in section "Preparatory tasks". The procedure is described under "Procedure" and the subsequent tasks are described under "Subsequent tasks".

**Troubleshooting / list of frequently asked questions:** – Clean the framework surfaces from grease and polishing residue. – Apply Sebond Grip only in thin layers. – Do not use intermediary polymerization units for polymerization.

Error	Cause	Remedy
does not cure	applied too thickly <p>insufficient polymerization</p>	apply only in thin layers <p>observe polymerization times</p> <p>– check lamp/replace lamp if necessary</p> <p>– extending the polymerization time will not effectuate a greater polymerization depth</p> <p>– do not use an intermediary polymerization lam</p>
material is too hard inside the syringe	syringe was stored at a temperature under 10° C (50° F)	store the material at room temperature or let the material reach room temperature in time before its use
greasy surface	insufficient polymerization	observe polymerization times <p>check unit/ regular maintenance</p>

**Working range of the material:** 1 - 3 min., depending on light conditions
Storage: Storage temperature 10 - 25° C / 50 - 77° F. Tightly close syringe immediately after each use.

**Shelf-life:** The maximum shelf-life is printed on the label of each packaging. Do not use after the expiration date.

**Side-effects:** With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (e. g. allergies) or local discomfort, however, cannot be ruled out completely. All serious incidents arising in connection with the use of this product must be reported to the manufacturer indicated below and to the relevant competent authority.

**Guarantee:** Our processing instructions, regardless of whether they are issued verbally, in writing or in the form of a practical demonstration, are based upon our own experience and can therefore be regarded solely as guidelines. Our products are subject to continuous development. We therefore reserve the right to make possible changes.

**Contra-indications / interactions:** If a patient has known hyper-sensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. In such cases, the composition of the medical device we supply is available upon request.

**Note:** Please supply the dentist with the above information, if this medical device is used to produce a special model. Please also consider the safety data sheets.

**Note on disposal:** Residual quantities and packaging material must be disposed of in accordance with local and/or statutory regulations.

**CE 0297**

**SCHÜTZ DENTAL**  
Micerium Group

Schütz Dental GmbH  
Dieselstr. 5 - 6 · 61191 Rosbach/Germany  
Telefon: +49 (0) 6003 814-0 · Telefax: +49 (0) 6003 814-906  
info@schuetz-dental.de · www.schuetz-dental.de