# Sebond Universal

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If there is anything in this instruction for use that you do not understand, please contact our customer service department before using the product.

As the manufacturer of this medical device, we inform our users and patients that all serious events occurring in connection with it must be reported to us (the manufacturers) as well as the relevant authorities in the Member State where the user and/or patient is resident.

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**Description / indication**

Sebond Universal is a one-component universal primer for the adhesive bond between a luting composite and indirect restorations made of lithium disilicate, glass ceramics, zirconium dioxide, high-precious alloys, non-precious alloys, titanium, Zirconia reinforced composite (PMMA-based), PEEK, PEKK.

**Contraindication**

The use of Sebond Universal is contraindicated  
 • if the patient has a known allergy against one of the product's components  
 • for applications outside of the indications described in this instruction leaflet.

**Composition**

Sebond Universal is an alcohol-based solution of silane methacrylate, phosphoric methacrylate and sulfide methacrylate.

**Application tips:**

- Close bottle immediately after each use.
- In case of a contamination of the restoration with saliva, blood or water, the complete treatment of the restoration must be repeated.
- Sebond Universal must not be taken out of the bottle directly with the brush, as this could lead to a contamination of the material inside the bottle.

**Instructions for use**

**Removal of the temporary restoration**

Remove the temporary restoration, clean the cavity or stump, remove residue of the temporary cement, flush well with water and subsequently, dry with oil-free compressed air.

**Trying in the restoration**

- Try in the restoration, check fit and contact points.
- Check the occlusion only very carefully in case of very brittle ceramic objects, as an object that has not yet been fixed by luting could crack or break.
- If necessary, corrections can be carried out with a fine diamond tool at medium speed and with only slight pressure. All surfaces that were ground must be polished again.
- Use alcohol (e.g. isopropanol) to thoroughly clean the objects to be fixed by luting and then dry them with oil-free compressed air.

**Note:** Adhere to the instructions of the framework material manufacturer.  
 • After drying the object, avoid moisture and contamination.  
 • For an adhesive fixation procedure, it is imperative to keep the working area dry, preferably with help of cotton rolls and gland patches. We recommend to keep the working area absolutely dry by use of a kofferdam.

**A) Fixation of indirect restorations**

We recommend the use of Capo Etch/Bond. It is also possible to use other commercially available dentin and enamel bonding agents.

**A1) Pre-treatment of the natural tooth**

- Apply Capo Etch to the preparation surface, etch enamel for 30 sec. and/or dentin for 15 sec.
- Afterwards, rinse off etching gel thoroughly with water.
- Wet the surface with a mixture of 1:1 CapoBond and CapoBond Activator.
- We recommend to apply the adhesive with a microbrush.
- Blow off water and solvent residue with a soft air stream.
- Apply a second coat of the mixture 1:1 CapoBond and CapoBond Activator and blow off again.

**A2) Pre-treatment of the restoration**

Generally follow the instructions of the manufacturer of the framework material.

Preparation of	Procedure
<ul style="list-style-type: none"> <li>• Metals</li> <li>• Zirconium dioxide</li> <li>• Aluminium oxide</li> </ul>	<ul style="list-style-type: none"> <li>• Sandblast the interior surfaces of the restoration with 30 - 50 µm aluminum oxide at max. 2 bar</li> <li>• Clean the restoration with alcohol and dry.</li> <li>• Evenly apply Sebond with a thin brush and let it sit for ca. 1 min. afterwards.</li> </ul> <p><b>Note:</b> In order to ensure an optimal bond, do not clean the surfaces with phosphoric acid!!!</p>
<ul style="list-style-type: none"> <li>• Lithium disilicate glass ceramics e.g. IPS e.max Press, IPS e.max CAD</li> </ul>	<ul style="list-style-type: none"> <li>• Etch with 5% hydrofluoric acid for 20 sec. or in accordance with the instructions of the framework material's manufacturer.</li> <li>• Thoroughly flush the restoration with water and dry with oil-free compressed air.</li> <li>• Evenly apply Sebond Universal with a thin brush and let it sit for ca. 1 min.</li> </ul>
<ul style="list-style-type: none"> <li>• PMMA materials, e.g. Tizian Zirconia Reinforced Composite (Schütz Dental)</li> <li>• PEEK</li> <li>• PEKK</li> </ul>	<ul style="list-style-type: none"> <li>• Sandblast the interior surfaces of the restoration with 30 - 50 µm aluminum oxide at max. 2 bar</li> <li>• Clean with a suitable cleaning agent, e.g. a disinfectant based on quaternary sal ammoniac. Afterwards, dry the objects.</li> <li>• Evenly apply Sebond Universal with a thin brush and let it sit for ca. 1 min. afterwards.</li> </ul>

**A3) Application of Sebond Universal**

Liberally apply Sebond Universal with a brush and let it sit for ca. 1 min. If necessary, apply a second layer.

**A4) Application of a luting composite**

Apply the luting composite (e.g. Alphaslink Cem) according to its instructions and insert the restoration.

**Contra-indications / interactions**

In case of irritations or if a patient has known allergies against or hyper-sensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. The dentist should consider known interactions and crossreactions of the product with other materials already in the patient's mouth before using the product. Unpolymerized resin may lead to allergic skin reactions. Therefore, the user should consider protective measures.

**Side-effects**

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (e.g. allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us. In order to avoid a possible pulp reaction, always prepare a liner in cases with exposed dentin (e.g. a compound containing calcium hydroxide).

**Storage**

- Sebond Universal must be stored at 10 - 25 °C (50 °F - 77 °F). After using the primer, close the bottle again immediately to keep the volatile components from evaporizing.

**Shelf life**

- The maximum shelf life is printed on the label of each primary packaging.
- Do not use Sebond Universal after the expiry date has passed.

**Note**

- Please consider the material safety data sheets that are available for this product.

**Only for use in dentistry!**