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If there is anything in this instruction for use that you do not understand, please contact our customer service department before using the product.

As the manufacturer of this medical device, we inform our users and patients that all serious events occurring in connection with it must be reported to us (the manufacturers) as well as the relevant authorities in the Member State where the user and/or patient is resident.

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# Instructions

## Intended Use

The millable CAD/CAM materials are used for the fabrication of dental prosthetic restorations with the aid of dental milling systems. The millable CAD/CAM materials are intended for use by the dental technician in a dental laboratory. Depending on the material used, the product can be used to create temporary or permanent dental restorations.

## General product description

The Tizian Blanks zirconia-reinforced composite are milling blanks made of zirconia-reinforced composite (PMMA-based, 20 % inorganic filler). They are processed in CAD/CAM procedures and can be used monolithically or veneered with composite. The blanks are available in various tooth shades.

## Composition

Mixed polymer based on methacrylate with 20 % inorganic filler as well as color pigments, initiators and stabilizers.

## Indications:

### 1. Long-term temporaries (for up to 2 years)

- a) Up to 16-unit anterior and posterior crowns and bridges, fully anatomic or reduced.
- b) Position no more than 2 bridge units in the posterior area.
- c) If the bridge is a free-end construction, add max. one posterior attachment.

### 2. Smaller permanent restorations\*

- a) individual crowns
- b) inlays
- c) onlays
- d) veneers
- e) bridges with a maximum of 3 units

\*It is highly recommended to veneer such restorations with an incisal material.

If the restoration is not veneered with a composite material, a high lustre polish (smooth surface) is indispensable.

## Contraindications

If a patient has hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. In such cases, we will supply the composition of our medical device upon request.

## Composition

> 78 % PMMA, 20 % inorganic filler, < 2 % inorganic pigments

## Safety note

During processing of „Tizian Blanks Zirconia Reinforced Composite“, dusts are released which can damage the lungs, and irritate the skin and the eyes. Therefore, please only process the material while running an adequate extractor system.

Please adhere to the following safety rules:

- Process the blanks only underneath a local extractor system.
- Ensure that product dust will not come in touch with the eyes or the mucous membranes.
- Do not inhale the dust.
- Wear a face mask and protective goggles.

## Interactions

The dentist should consider known interactions and crossreactions of the product with other materials already in the patient's mouth before using the product.

## Design notes / framework dimensioning

- Minimum wall thickness:
  - 0.6 mm (cervical)

- 1.0 mm (occlusal)
- > 0.2 mm (crown margin)
- Diameter of bridge units: max. 10 mm
- Profile of connectors:
  - > 10 mm<sup>2</sup> (anterior)
  - > 15 mm<sup>2</sup> (posterior)

Important: To increase the stability of the construction, choose the greatest possible height for the connectors.

### Milling procedures

- The optimal ambient temperature for milling this material is 20 °C – 25 °C (68 °F – 77 °F).
- Please adhere to the instructions of the milling system. We recommend to use **a solid 5-axes milling system** to ensure an excellent fit of the restorations.
- Tizian Blanks Zirconia Reinforced Composite are suited to dry milling.
- All milling **tools** must be suited to processing composites. We recommend coated radius cutters which are conventionally used for composites.
- The choice of **milling parameters** strongly depends on the milling system, on the CAM-software and on the milling geometries used. If you are not certain which parameters to use, please use your parameters for composites.

### Detaching the restoration

Carefully detach the finished restoration from the blank with a jig cutter.

### Finishing / polishing

- Use clean cross-cut hardmetal mills to finish the restorations.
- Use a low number of revolutions and a minimum of pressure to avoid thermal damage. Ensure sufficient cooling.
- Polish with conventional pastes for polishing composites (e. g. dialog Vario

polish), goat bristle brushes and cotton polishers.

### Important:

- Do not overheat the material when using mills, steam cleaners or similar.
- The wall thickness of the frameworks must not go below a limit of 0.4 mm.

### Cleaning the frameworks

- Alcohol in high concentrations will generally cause cracks on PMMA, the so-called crazing effect. Do not use cleaners or disinfectants which contain alcohol.
- Do not overheat the surface of the restoration while steam-cleaning it.
- Alternative possible cleaning methods are mechanical procedures, using an ultrasonic unit or alcoholfree disinfectants with quaternary ammonium compounds.

### Veneering

After suitable surface conditioning, Tizian Blanks Zirconia Reinforced Composite can be veneered with all conventional veneering composites:

### Recommended procedure

1. Condition the surfaces by sandblasting with aluminium oxide (max. 110 µm and 2 bar, obtuse angle). Ensure that the surfaces are clean, dry and grease-free, before you continue processing the material.
2. Apply a suitable bonding agent, e. g. dialog Bonding Fluid (Schütz Dental).
3. Veneer conventionally, e. g. with dialog Vario Composite Paste and/or Flow. Observe the instructions of the composite's manufacturer.

Important: When preparing a vestibular veneer, it is necessary to apply a protection on the masticatory edge.

**Tip:**

If you are using the Schütz Dental veneering system, dialog Occlusal masses can be applied directly to the framework, without a foundation of paste or flow material. To obtain an even bond in such cases, we recommend the use of Sebond Grip.

When using a different veneering system, the incisal mass should not be applied directly to the framework. This could lead to chipping, depending on the material used.

**Fixation**

The material must only be fixed by adhesion with a conventional temporary or permanent fixation composite. Please adhere to the instructions of the fixation composite's manufacturer.

**Notes**

- Do not use cements of any kind.
- Please pass on all the above information to the attending physician if you are processing this medical device for a custom-made product. Observe the relevant safety data sheets during processing.

**Storage**

Store in a dry place at temperatures between 10°C and 25°C (50°F and 77°F). Keep away from direct sunlight and moisture.

**Shelf-life**

The maximum shelf-life is printed on the label of the respective packing. Do not use after the shelf-life has expired.

**Side effects**

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (e. g. allergies) or local discomfort, however, cannot be ruled out completely.

Any serious incidents arising in connection with the use of this product must be reported to the manufacturer indicated below and to the relevant competent authority.

**Note on disposal**

Residual quantities and packaging material must be disposed of in accordance with local and/or statutory regulations.

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