Tizian occlusal

CE 0297

SCHÜTZ  **DENTAL**
Micerium Group

Schütz Dental GmbH · Dieselstr. 5-6 · 61191 Rosbach/Germany
Tel.: +49 (0) 6003 814-0 · Fax: +49 (0) 6003 814-906
info@schuetz-dental.de · www.schuetz-dental.de

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If there is anything in this instruction for use that you do not understand, please contact our customer service department before using the product.

As the manufacturer of this medical device, we inform our users and patients that all serious events occurring in connection with it must be reported to us (the manufacturers) as well as the relevant authorities in the Member State where the user and/or patient is resident.

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Processing information

Description / Indication

Tizian occlusal is a translucent zirconium oxide material used for the production of full zirconium bridges. The material is specifically adapted for colouring techniques without ceramic veneering. Tizian occlusal can be used particularly where space is limited, in implant areas or treatments involving gingival construction. Tizian Occlusal can be fully anatomically formed as well as classically veneered.

Safety Precautions

Tizian occlusal is delivered in the form of white-fired blanks. The blanks possess only a limited strength and display a residual porosity. Therefore, careful handling is essential. Please check the condition of the packaging and product directly after receipt of the delivery.

Attention: The pre-sintered structure must never come into contact with water or steam.

General information

The information provided in this instruction manual must be passed on to any person using the products mentioned therein. The products must only be used by qualified personnel. The user is obliged to use the products in accordance with the present instruction manual and with appropriate hygiene measures and to verify on his/her own responsibility whether the products are suitable for the individual patient situation. The user will be held fully responsible for the appropriate and correct use of the products. The manufacturer assumes no liability for incorrect results in form of direct or indirect damages or any other damages that occur from the use and/or the processing of the products. Any claim for damages (including punitive damages), is limited to the commercial value of the products. Independently of this, the user is obliged to report all serious incidents that occur in connection with

the products to the competent authority and to the manufacturer.

Contraindications

Should intolerance to the products be suspected, use them only after an allergy screening test has been performed with issue of medical statement of non-allergy conditions.

Processing of Tizian occlusal

1. Preparation guidelines

Preparations suitable when using zirconium oxide are the groove or shoulder preparations. Unsuitable preparations, e.g. gutter or tangential preparations can lead to failure. The horizontal angle of preparation should be 4–6° and the vertical angle 3–5°. All margins should be rounded off and smoothed, avoid any sharp angles. Be sure of an optimum insertion between the abutments. The occlusal ablation should be at least 1.5 mm in the posterior tooth area and 2 mm in the anterior tooth (incisor) area.

2. Form-grinding

To produce a dental prosthesis from Tizian occlusal only those machines, equipment and tools should be used which are approved for working white-fired and presintered zirconium oxide material.

During processing, the white body must not come into contact with the cooling fluid. Please adhere to the instructions in your miller operating manual.

3. Framework dimensioning

The wall thickness before sintering must be at least 0.7 mm. After densification sintering the wall thickness of the crown must not fall below 0.4 mm. The wall thickness at the edge of the crown should not be less than 0.2 mm.

Dimensioning of connectors

The products of class 4 according to ISO 6872 are intended for bridges with a max. span of up to 3 units. When manufacturing bridges, the dimensions of the connectors need to be taken into account. It should be 3 x 4 mm in the anterior region and 4 x 4 mm in the posterior region.

The geometry of the connector is of paramount importance for the crack resistance of frameworks produced from zirconium oxide. Therefore, the selected height of the connector must be as large as possible.

In the posterior tooth area, a maximum of 2 pontics between the abutments is permitted. For cantilever bridges a maximum of 1 end pontic is allowed. In general: the more pontics required the more stable the framework needs to be.

In order to avoid so called „chipping“, anatomical modelling of the framework should be carried out so that an even veneer ceramic layer of 1–2 mm thickness can be applied.

4. Colouring

All the usual commercially-available colouring fluids for zirconium oxide can be used for colouring the frameworks, following the manufacturer's instructions. Here it is particularly important to adhere to the drying times. The best colouring results are achieved when the milled work-piece is not dipped into the colouring solution, but rather hand-painted using a brush. This procedure is comparable with the colouring technique before the final firing. In order to prevent discolouration, the brush must not have metal fittings, instead should be simply from plastic or wood. Please observe the instructions provided by the manufacturer of the colourant.

5. Sintering in a high temperature oven

Before the densification sintering, the milling dust needs to be removed from the frameworks. The 25% larger milled

framework will shrink to the required size during the sintering process. Thereby the components are conferred their final properties. A suitable furnace is essential for the sintering process, one that is approved for the sintering of white-fired zirconium oxide components.

Sintering programme:

Start: room temperature

Heating rate: 8° C/min.

End temperature: 1600° C

Holding time at end temperature: 120 min.

Cooling rate: 8° C/min.

Opening of the oven at 300° C

It is imperative that large bridges are supported. Only in this way can a distortion due to sintering be avoided. Zirconium oxide is a poor heat conductor. As a result, the zirconium oxide in solid pontics is significantly harder to bake than that in thin-walled caps. Therefore, we recommend the heating and cooling rates to be reduced to 5° C/min. for large pieces and large volumes of work.

Chart of sintering parameters – heating rates

Weight (incl. sintering stabiliser / dental)	Heating rate	Holding time
0 – 2 g	max. 14° C/min.	1,5 h
2 – 4 g	max. 8° C/min.	2 h
> 4 g	max. 5° C/min.	2 h

6. Finishing the frameworks

The finishing touches to the unsintered work are performed by hand. Zirconium oxide stones, sinter diamonds, rubber polishers and hard metal millers are suitable.

The dense sintered frameworks are preferably processed non-mechanically.

If this cannot be avoided, then only water-cooled, diamond inset tools in perfect condition should be used, otherwise the structure of the pieces could be damaged. Too intense localised heating must be

avoided (e. g. by steam blasting, sand blasting or polishing to a shine).

High pressure should not be applied during processing. The thinner the wall, the greater the care required. Care must be taken to ensure the wall thickness does not fall below 0.4 mm. If possible the interdental connection points should not be ground. Grooves or grinding must on no account be made at the basal connection points. Avoid any sharp edges. The use of cutting discs should be avoided. Sand blasting is not recommended. Should this be necessary, the framework can be blasted after conditioning with aluminium oxide at a max. of 120 µm at < 2 bar pressure. Afterwards the framework should be washed with water.

Should the sintered framework require a lot of grinding, a repair bake must be carried out: heating rate 8° C/min., end temperature 1055° C, holding time 5 min., natural cooling.

7. Facing

The facing must be exclusively realised using a facing ceramic suitable for zirconium oxide, e. g. Tizian. The manufacturer's handling instructions must be observed.

In order to reach an optimal adhesion, we recommend performing a wash bake with Dentin (100° C higher than the normal bake temperature).

Zirconium oxide is a poor heat conductor. As a result, the facing ceramic in solid pontics is significantly harder to bake than on crowns. Therefore there is a special baking procedure.

If the bridge to be faced has no solid pontics, we recommend the usual heat increase rate of 55° C. However, should the volume of the pontics rise, this requires an increase of 2 mins. in the holding time at the end temperature, as compensation for the different heat con-

ductivity in the crown and the pontic. Here it is recommended that the temperature be reduced by around 10° C, to avoid excessive ceramic melting. In general, during the use of zirconium oxide layered ceramic, a two-minute holding time under vacuum at the end temperature should be adhered to.

Attention: The more massive the framework is, the lower the heating rate that should be selected. Furthermore, slow cooling-off should be used to avoid the formation of cracks. A minimum holding time of 2–3 minutes must be adhered to. Temperature shocks while baking (especially with massive frameworks) must be avoided. The framework should only be taken out of the furnace when the temperature has dropped below 200° C. Warm work must not be placed on a cold surface (such as a table top) as this can lead to cracking.

For the veneering of bridges of various elements you must observe the following parameters during the heating and the cooling cycles:

Weight (incl. sintering stabiliser / dental)	Heating rate	Holding time
0 – 1 g	per dental unit	55° C/min.
1 – 2 g	per dental unit	45° C/min.
2 – 3 g	per dental unit	35° C/min.
über 3 g	per dental unit	25° C/min.

8. Fitting

The work can be fitted conventionally with phosphate or glass ionomere cement and adhesive attachment is also possible.

For the incorporation of the manufactured restorations, only materials may be used that are approved and tested for zirconium oxide. The manufacturer's handling instructions must be observed.

Important information

When working with partly and densely sintered zirconium oxide, dust can be released,

which can be absorbed through the alveoli, can damage the lungs, as well as lead to irritation of the eyes and the skin. Therefore it must only be processed when the ex-hauster is in good working order. The technicians should wear eye protection and a dust mask.

Storage

The blanks must be stored under dry conditions at room temperature. The white bodies must not be exposed to impacts or vibration. In general, contamination is to be avoided.

Do not expose to heavy shocks or vibrations. Blocks must not come into contact with liquids. Avoid any exposure to humidity. Avoid the formation and deposition of dust.

Disposal

The material can be disposed of in the household waste.

Liability

Our advice for use, either in written, verbal or as practical guidance is based on our own experience and therefore can only be seen as guidelines.

Side effects

Unwanted side effects of this medical product are extremely rare when the product has been handled and used appropriately.

Immune reactions (e. g. allergies) or local discomfort cannot be ruled out.

Should you be made aware of any unwanted side effects – even in case of doubt – we request that we are informed.

Contraindications / interactions

If a patient is hypersensitive to one or more of the components, this product may not be used or only be used under strict supervision of the treating surgeon / dentist. Known cross-reactions or interactions of this medical device with other materials already present in the mouth must be taken into account by the treating dentist when used.

Important information

Please forward all of the above information to the treating dentist if you are processing this medical device for a special design. When working, please adhere to the specific product safety sheets.

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