dialog Vario Chroma-Flow

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Light-curing stains on composite basis

Enthält Aliphatisches Urethandimethacrylat, 1,4-Butandioldimethacrylat, Diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide

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Contains di-urethane dimethacrylate, tetramethylene dimethacrylate, diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide

Warning. May cause an allergic skin reaction. Toxic to aquatic life with long lasting effects. Avoid breathing vapours / spray. Avoid release to the environment. Wear protective gloves. If skin irritation or rash occurs: Get medical advice / attention.

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If there is anything in this instruction for use that you do not understand, please contact our customer service department before using the product.

As the manufacturer of this medical device, we inform our users and patients that all serious events occurring in connection with it must be reported to us (the manufacturers) as well as the relevant authorities in the Member State where the user and/or patient is resident.

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Instructions for use

Range of application

Characterization of composite facings on crowns and bridges, individualization of tailor-made resin teeth as well as characterization of denture base resins in partial and full dentures, individualization of fixed restorations inside the patient's mouth. Inlays, Onlays, Veneers and fillings.

Available in ten different shades.

Instructions

The stains can be applied directly from the syringe. Alternatively, the stains can be put on a palette for mixing which can be covered and can be applied with a thin brush. All shades are intermixable.

Consider the instructions for the composite.

Apply dialog Vario Chroma-Flow in layers no thicker than 0.3 mm. We recommend the placement of a rubber dam when using the material inside the patient's mouth.

Polymerization

dialog Vario Chroma-Flow can be polymerized in all conventional light-curing systems.

Before the final polymerization, the complete polymerization site may be covered with expers-Gel. This avoids the formation of a new dispersion layer and facilitates the finalization of the object.

Use a light-curing system with an emission spectrum of 310 - 500 nm. Perfect physical properties can only be obtained with immaculate lamps. It is important to test the light intensity regularly, in line with the manufacturer's instructions. When using a handheld polymerization lamp, it is important to ensure that the complete coat of material is exposed to the light.

Finishing

All silicone polishers (square-edge or knife-edge wheels, cylinder) as well as hard metal and diamond mills are suitable for finishing the material.

The veneer is polished with handpiece-mounted goat-hair brushes, woolen wheels and a polishing paste.

Careful conditioning of the surface is essential if optimal results are to be achieved, thus virtually precluding deposits (nicotine, caffeine etc.) which could cause discoloration.

Working range of the material

1 - 3 minutes, depending on the light conditions.

Storage

To retain the material's stability (shelf-life), store at temperatures of 10 - 25 °C or 50 - 77 °F. Close syringe directly after use.

Shelf-life

The maximum shelf-life is imprinted on each syringe. Do not use after shelf-life has expired.

Side-effects

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us.

Contra-indications / interactions

If a patient has known allergies against or hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. The dentist should consider known interactions and cross reactions of the product with other materials already in the patient's mouth before using the product.

Note

Please supply the dentist with the above information, if this medical device is used to produce a special model. Please also consider the safety data sheets.

Interactions with other materials

As phenolic substances (such as eugenol) inhibit polymerization, do not use materials containing these substances.

Disposal

Must be disposed of in a special waste deposit in accordance to government regulations.

Composition

Filler, anorganic (50 weight % respectively 36% by volume)

Glass filler (medium particle size: 0.7 µm)

Pyrogenic silica (medium particle size: 0.01 µm)

Monomers (49 weight %)

Di-urethane dimethacrylate

Tetramethylene dimethacrylate

Bis-GMA

Additives (1 weight %)

Initiators, stabilizers, pigments

Polymerization times

Spektra LED	1 min.	Spektramat	5 min.
Spektra 2000 fast	3 min.	Labolight LV- II / III	5 min.
Spektra 2000	5 min.	Solidilite EX	5 min.
Hera Flash / HiLite	3 min.	Bluephase	20 sec.